


# Comfort Plus™ Premium Saliva Ejectors


## Make your patient's dental experience more comfortable.

Comfort Plus™ Premium Saliva Ejectors are engineered with long, rounded edges, smooth surfaces, and strategically placed suction ports with curved openings to protect delicate mucosal tissue and optimize fluid and light debris removal.

- **Increases Patient Comfort:** Suction ports positioned along the side rather than wrapping around the end of the tip to minimize trauma to mucosal tissue.
- **Prevents Blocked Suction:** Engineered to optimize suction for uninterrupted fluid and light debris removal.
- **Sturdy and Robust:** Holds shape when formed and placed in mouth.



Perfectly placed suction ports take in fluid, not tissue.



### Comfort Plus™ Premium Saliva Ejectors

|                |                                           |         |
|----------------|-------------------------------------------|---------|
| <b>ZWWCP</b>   | White with White Tip                      | 100/Bag |
| <b>ZCWCP</b>   | Clear with White Tip                      | 100/Bag |
| <b>ZWGPCM</b>  | White with Green Tip<br>Mint Scented      | 100/Bag |
| <b>ZWPCPBG</b> | White with Pink Tip<br>Bubble Gum Scented | 100/Bag |

Call your preferred distributor to order, or learn more at [HuFriedyGroup.com](https://www.HuFriedyGroup.com)

Not made with natural rubber latex.

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions. [HuFriedyGroup.com](https://www.HuFriedyGroup.com)  
©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-842 / 2607