



A SHARPER EDGE FOR BETTER CARE

| EverEdge® Technology

How the best perform



EVEREDGE® TECHNOLOGY

Hu-Friedy's EverEdge Technology is unlike anything you've experienced before in a scaler. We've combined state-of-the-art technology in metallurgy, heat treatment and cryogenics to create a superior stainless steel that keeps your scalers and curettes sharper longer. That means less frequent sharpening, less hand fatigue, and greater comfort throughout the day.

“As an educator, I always recommend Hu-Friedy's EverEdge instruments to my dental hygiene students. When well maintained, they can be sharpened so quickly and easily. A real advantage of the fine sharp cutting edge is **less hand fatigue for the hygienist and better comfort for the patient**, so both the clinician and patient can feel the benefits of EverEdge.”

— Jill Nield Gehrig, RDH, MA
Dean Emeritus, Division of Allied Health

STAYS SHARPER LONGER

The revolutionary stainless steel alloy is super-durable and stays sharper longer so the clinician can spend less time sharpening.

PROPRIETARY PROCESS

Our proprietary processing ensures that the superior edge retention and wear characteristics of EverEdge Technology will last the entire life of the instrument. Since it's not a superficial coating, they can be resharpened like any other scaler to extend the useful life.

FINELY-HONED FINISH

EverEdge Scalers feature a hand-crafted finish from Hu-Friedy's professional artisans that keeps edges effectively sharp and allows for efficient scaling.

DIAMOND KNURL


The unique diamond knurl pattern on the large diameter handle provides a confident and comfortable grasp. The #9 EverEdge handle was also designed to be the optimal weight for reduced hand fatigue.

COMFORT ZONE CONE WITH SIGNATURE SERIES GRIP

The cushioned grip provides a smooth transition and increased comfort. With distinctive Ocean Blue grips, EverEdge scalers can be easily identified in a set up or customized for your office's needs.

NOW AVAILABLE IN RESIN 8 COLORS

These lightweight resin handles use a color-coded system to easily identify scalers and curettes that correspond to anatomical areas of the mouth.



SCALER CHECK-UP

Hu-Friedy uses the highest quality materials and processes to create EverEdge scalers' long lasting cutting edges but they still need to be regularly maintained and then replaced after a certain amount of time. The useful life of a scaler depends on several factors including: frequency of use, frequency of sharpening, sterilization method and clinical applications. Hu-Friedy recommends that you regularly audit your scalers to determine when they should be replaced.

AUDIT CHECKLIST:

- ☐ Retipped instruments
- ☐ Dull/chipped cutting edge or point
- ☐ Broken or bent tips
- ☐ Thin/worn-out blades
- ☐ Pointed curettes
- ☐ Corrosion

For more information about
Hu-Friedy's Total Solutions visit


HU-FRIEDY.COM

TOP 10 SCALERS

1	H6/H7	SH6/79
2	204S	S204S9
3	Nevi 4	SCNEVI49
4	Gracey 13/14	SG13/149
5	Gracey 11/12	SG11/129
6	Columbia 13/14	SC13/149
7	Barnhart 5/6	SBH5/69
8	Gracey 1/2	SG1/29
9	Nevi 2	SCNEVI29
10	H5/33 Jacquette	SH5/339

DID YOU KNOW?

Curettes are designed with a rounded toe for patient comfort during subgingival scaling. Sharpening can turn this "toe" into a point and cause patient discomfort when used subgingivally.



PERIODONTAL



RESTORATIVE



ORTHODONTIC




INSTRUMENT MANAGEMENT



ENDODONTIC



SURGICAL



ULTRASONIC SCALING



Hu-Friedy Mfg. Co., LLC: 2501 Sanders Rd., Northbrook, IL 60062 | USA
Telephone: 1-800-Hu-Friedy (1-773-975-6100) | Website: Hu-Friedy.com

©2026Hu-Friedy Mfg. Co., LLC. All rights reserved. HF-479/0914

