


Black Line Endo Collection

Designed for Premium Performance




Instrumentation and
Procedural Products



Endodontic Instrumentation

Our comprehensive line of endodontic instruments and accessories give you versatile options to meet all of your clinical needs. If you are looking for exceptional instrument design, quality and reliability, and unmatched customer service, choose the HuFriedyGroup Endodontic Instruments.



Learn more about
Endodontic Instrumentation

THE CHALLENGE

The sensitive nature of endodontic procedures necessitates instruments with exceptional response, integrity and durability to accommodate a wide variety of endodontic procedures—from routine root canals to complex microsurgery.

THE SOLUTION


Our Black Line Endodontic Instruments are individually handcrafted to precise specifications from quality materials to ensure superior performance, longevity and reliability. Instrument design combined with proprietary materials and incomparable workmanship result in instruments with unique balance, working leverage and control. They function as a seamless extension of the practitioner's skill.

THE RESULT

High-quality HuFriedyGroup Instrumentation provides unsurpassed clinical efficiency, so practitioners can deliver superior outcomes and focus on the patient. The result is outstanding patient care.


Black Line for Endo

Featuring a performance engineered coating to enhance lubricity, reduce light reflection under harsh lighting, and optimize edge retention, the Black Line Endo Instruments are designed to provide exceptional clinical outcomes and improved efficiency during endodontic procedures.



Explorers

Our endodontic explorers are made from a special blend of alloy which provides the flexibility, strength, tactile sensitivity and resilience clinicians need when exploring and assessing the treatment site. They are available in several double-ended designs and handles.




DG16
EXDG16X

DG16/17
EXDG16/17X

DG16/23
EXDG16/23X

Excavators

Our endodontic excavators feature extended-length shanks and are meticulously handcrafted from the highest grade stainless steel. They help the clinician access the pulp chamber for easy and efficient removal of dentin. The result is a superior clinical outcome and improved efficiency.



Spreaders

We manufacture a wide variety of spreaders, all handcrafted from the finest stainless steel to precise specifications. Our spreaders help the clinician gain access into the most difficult canals for optimal compaction of filling material during lateral condensation.



D11TS Thinnest (21 mm)
RCSD11TSX



MA57 (27 mm)
RCSMA57X



GP1 (20 mm)
RCSGP1X




D11T Thinner (22 mm)
RCSD11TX



D11 Standard (23 mm)
RCSD11X

Pluggers

Our pluggers aid clinicians in performing vertical condensation as they are handcrafted to exact specifications using the finest materials.






Berutti Pasqualini Plugger
0.5 mm & 0.7 mm
MBP-PLG1



Berutti Pasqualini Plugger
0.5 mm
MBP-PLG2



Berutti Pasqualini Plugger
0.7 mm
MBP-PLG3




Black Line Rubber Dam Clamps

Our rubber dam clamps are made in the USA from proprietary Satin Steel™ Alloy for corrosion resistance, superior strength, and a fine balance between flexibility and memory. The black matte finish absorbs and dissipates light instead of reflecting it back, so you can experience a comfortable field of vision that is easier on the eyes—to help you be the Best in Practice.

SKU	DESCRIPTION
RDCM106LX	Rubber Dam Clamp #106L, Black Line
RDCM106RX	Rubber Dam Clamp #106R, Black Line
RDCMB4X	Rubber Dam Clamp #B4, Black Line
RDCM12AX	Rubber Dam Clamp #12A, Black Line
RDCM13AX	Rubber Dam Clamp #13A, Black Line
RDCM14AX	Rubber Dam Clamp #14A, Black Line
RDCM2AX	Rubber Dam Clamp #2A, Black Line
RDCM14X	Rubber Dam Clamp #14, Black Line
RDCM9X	Rubber Dam Clamp #9, Black Line
RDCM2X	Rubber Dam Clamp #2, Black Line
RDCM00X	Rubber Dam Clamp #00, Black Line
RDCM7X	Rubber Dam Clamp #7, Black Line
RDCM212X	Rubber Dam Clamp #212, Black Line
RDCM8AX	Rubber Dam Clamp #8A, Black Line
RDCMW8AX	Rubber Dam Clamp #W8A, Black Line
RDCM1X	Rubber Dam Clamp #1, Black Line
RDCM0X	Rubber Dam Clamp #0, Black Line
RDCM8X	Rubber Dam Clamp #8, Black Line
RDCM27NX	Rubber Dam Clamp #27N, Black Line
RDCMW3X	Rubber Dam Clamp #W3, Black Line
RDCM3X	Rubber Dam Clamp #3, Black Line
RDCM27X	Rubber Dam Clamp #27, Black Line
RDCMW14AX	Rubber Dam Clamp #W14A, Black Line
RDCM26NX	Rubber Dam Clamp #26N, Black Line
RDCM4X	Rubber Dam Clamp #4, Black Line
RDCM6X	Rubber Dam Clamp #6, Black Line
RDCMW2X	Rubber Dam Clamp #W2, Black Line

SKU	DESCRIPTION
RDCM2ASX	Rubber Dam Clamp #2AS, Black Line
RDCM56X	Rubber Dam Clamp #56, Black Line
RDCM201X	Rubber Dam Clamp #201, Black Line
RDCM18X	Rubber Dam Clamp #18, Black Line
RDCM7AX	Rubber Dam Clamp #7A, Black Line
RDCM56SX	Rubber Dam Clamp #56S, Black Line
RDCM205X	Rubber Dam Clamp #205, Black Line
RDCM5X	Rubber Dam Clamp #5, Black Line
RDCM8ADX	Rubber Dam Clamp #8AD, Black Line
RDCM9SX	Rubber Dam Clamp #9S, Black Line
RDCM212SAX	Rubber Dam Clamp #212SA, Black Line
RDCM1AX	Rubber Dam Clamp #1A, Black Line
RDCMW7X	Rubber Dam Clamp #W7, Black Line
RDCM209X	Rubber Dam Clamp #209, Black Line
RDCMW4X	Rubber Dam Clamp #W4, Black Line
RDCM26X	Rubber Dam Clamp #26, Black Line
RDCM24X	Rubber Dam Clamp #24, Black Line
RDCM25X	Rubber Dam Clamp #25, Black Line
RDCM29X	Rubber Dam Clamp #29, Black Line
RDCM10X	Rubber Dam Clamp #10, Black Line
RDCM28X	Rubber Dam Clamp #28, Black Line
RDCM11X	Rubber Dam Clamp #11, Black Line
RDCM31X	Rubber Dam Clamp #31, Black Line
RDCM30X	Rubber Dam Clamp #30, Black Line



RDCOB

Rubber Dam Clamp Organizing Board	RDCOB
IMS™ Parts Box	IMS-1271
IMS™ Parts Box	IMS-1273



RDCM106LX, RDCM106RX, and RDCMB4X



Find the perfect clamp for your procedure

For more information visit HuFriedyGroup.com/BlacklineRDC

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | HuFriedyGroup.com
All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.
©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. HF-656/1023