

ARTICULATING & OCCLUSION TESTING MATERIALS



- Combined use of articulating film and paper to mark high and low occlusion points provides the most accurate way to mark every point of articulation.
- Proprietary manufacturing process combined with stringent quality control guarantees a clear, positive marking without smudging of wet or dry fields.
- Color-coded packaging and product descriptions make identification easy.
- Store in dry, room temperature locations.



Articulating Paper 12 sheets/book - 12 books/box

TPXXT	XX Thin .00125"/32 microns	Blue
TPXT	Extra Thin .0015"/38 microns	Blue
TPT	Thin .0028"/71 microns	Blue
TPBR	Combo .0028"/71 microns	Blue/Red
TPTH	Thick .005"/127 microns	Blue
12 sheets/book - 6 books/box		
TPH	Horseshoe .0035"/89 microns	Blue/Red



Articulating Film 25 sheets/book - 12 books/box .0005"/13 microns

TPP	Precisely True Thin	Blue
-----	---------------------	------

Learn more at [HuFriedyGroup.com](https://www.HuFriedyGroup.com)