

HuFriedyGroup Key Opinion Leader

Dr. Robert Stanley D.D.S., PhD, MS

Extraction Kit



Focus on: **Implantology**



Country: **USA**



Website:
**stanleyinstitute.com/
doctors/dr-robert-stanley/**



Instagram: **@stanleyinstitute**



Facebook: **Facebook.com/
stanleyinstitute**



YouTube:
@stanleyinstitute

Dr. Robert Stanley D.D.S., PhD, MS

Dr. Stanley's professional career began as a Mechanical & Aerospace engineer. After receiving his BS, MS and PhD from North Carolina State University in engineering with a concentration in control theory, he enjoyed almost a decade with a telecommunications company in Research Triangle Park. To pursue his passion of dentistry, he returned to school and earned his dental degree from the University of North Carolina at Chapel Hill School of Dentistry, where he became Smile Engineer™.

As a dentist and an engineer, Dr. Robert Stanley has the unique ability to master and share his understanding of biomechanics as it applies to the daily implementation of dental implants. Dr. Stanley has extensive experience in implant surgery, as well as advanced bone and soft tissue grafting.

He is a Diplomate with the American Board of Oral Implantology, the International College of Oral Implantologists (ICOI), certified in Oral, IV and Conscious Sedation and is a member of the Academy of General Dentistry, the American Academy of Cosmetic Dentistry (AACD), the American Dental Association (ADA), the Raleigh-Wake Dental Society, the North Carolina Dental Society, the Dental Organization for Conscious Sedation (DOCS) and the International Association of Comprehensive Aesthetics.

This specialized surgical kit features Dr. Stanley's complete surgical extraction set-up. This instrumentation is designed for applications in simple and surgical atraumatic extractions and extractions with immediate implant guided placement. Additionally, this kit includes instrumentation for suturing.



Stanley Extraction Kit Components

Oral Surgery Signature Series™ Cassette, Yellow	IMOS5
Cone Socket Mirror Handle, Black	MH6B
5 HD Cone Socket Mirror Single Sided, Black	MIR5HDB
TU-17/23 Double Ended Explorer	EXTU17/236
2 Curved Kelly Hemostat	H2
University of Minnesota Cheek Retractor	CRM
9 Molt Periosteal Elevator, Black Line	P9X
5 Straight Round Scalpel Handle	10-130-05
87 Lucas Surgical Curette* (QTY: 2)	CL876
General Perma Sharp™ Tissue Pliers, 6"	TP5070
General 1x2 Perma Sharp™ Tissue Pliers, 6"	TP5071
13/14 Heidbrink Double Ended Root Tip Pick, Black Line	EHB13/14X
1X Miller-Colburn Bone File	BF1X6
9 Dean Scissors	S9
Straight Iris Perma Sharp™ Scissors	S5082
301 Apexo Elevator	E301
34S Straight Seldin Elevator	E34S
1L Seldin Elevator	E1L
Castroviejo Perma Sharp™ Needle Holder	NH5024R
CW Aspirating Anesthetic Syringe	SYRCW

Recommended Kit Additions (sold separately):

Bone Mill, Titanium	BMTI
Beyer Rongeurs	RBEY
#18 DE Ball Burnisher	BB18
#DG16 DE Endo Explorer, Satin Steel	EXDG166
Molt Mouth Gag, Adult Size	MGA

For more information, visit **HuFriedyGroup.com**

*Custom Serrated, serration service fee added

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | HuFriedyGroup.com

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.

©2024 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-693/2409

HuFriedyGroup
The Best In Practice