


# Disposable Prophy Angles

## Twist™ Plus Prophy Angles



## POINTS OF PERFORMANCE

- The Twist™ Plus Prophy Angle is the ONLY disposable that oscillates 120° in a smooth, continuous counter-rotational motion to eliminate the splatter and frictional heat associated with traditional 360° spinning prophy angles.
- Single-use design for infection control compliance.
- Round head design improves visibility inside the oral cavity and increases maneuverability and patient comfort.
- Superior cup design with external fluting: Enhanced inner web provides better pickup and retention of paste, while external fluting improves interproximal stain/biofilm removal.
- Recessed ridges in the grip makes a modified pen grasp more comfortable and secure while wearing gloves.
- Handpiece connection is on the underside, which reduces the potential for hand movements to push the angle off the handpiece.
- Tight seal of the cup to head reduces the risk of catching soft tissue in the corners of the mouth.

| SKU              | Description                    | Quantity |
|------------------|--------------------------------|----------|
| <b>TPLUSPASC</b> | Twist™ Plus DPA, Soft Cup      | 100/box  |
| <b>TPLUSPAFC</b> | Twist™ Plus DPA, Firm Cup      | 100/box  |
| <b>TPLUSPATB</b> | Twist™ Plus DPA, Tapered Brush | 100/box  |

# Disposable Propphy Angles

## Twist™ Propphy Angles




### POINTS OF PERFORMANCE

- Ergonomically designed angle that provides a 90° reciprocating “back and forth” rotating motion for continuous contact.
- Smooth back and forth motion reduces frictional heat, keeping the tooth surface cooler and reduces splatter.
- Soft web cup and ribbed flair edges, provides easy application and spreading of prophylactic paste.
- Single-use.
- Propphy angles not made with natural rubber latex.
- Made outside of the USA.

| SKU           | Description                         | Quantity |
|---------------|-------------------------------------|----------|
| <b>TPAFC</b>  | Twist™ Propphy Angle, Firm          | 100/box  |
| <b>TPASC</b>  | Twist™ Propphy Angle, Soft          | 100/box  |
| <b>TPATB</b>  | Twist™ Propphy Angle, Tapered Brush | 100/box  |
| <b>TPAFB</b>  | Twist™ Propphy Angle, Flat Brush    | 100/box  |
| <b>TPASCB</b> | Twist™ Propphy Angle, Soft          | 500/box  |

## Sparkle™ Propphy Angles



### POINTS OF PERFORMANCE

#### Created with hygienists' feedback:

- Patented beveled gear design runs vibration-free/ noise-free.
- Flared cup design improves subgingival access and reduces spatter, minimizing cross-contamination.
- Light colored angle with dark colored cup provides greater visibility.
- Single-use.

| SKU            | Description                           | Quantity |
|----------------|---------------------------------------|----------|
| <b>TPAFEZ</b>  | Sparkle™ EZ Contra Angle, Firm        | 100/box  |
| <b>TPASEZ</b>  | Sparkle™ EZ Contra Angle, Soft        | 100/box  |
| <b>TPASFC1</b> | Sparkle™ Propphy Angle, Firm          | 100/box  |
| <b>TPASSC1</b> | Sparkle™ Propphy Angle, Soft          | 100/box  |
| <b>TPATBC1</b> | Sparkle™ Propphy Angle, Tapered Brush | 100/box  |
| <b>TPASFC5</b> | Sparkle™ Propphy Angle, Firm          | 500/box  |
| <b>TPASSC5</b> | Sparkle™ Propphy Angle, Soft          | 500/box  |

For more information, visit [HuFriedyGroup.com](https://www.HuFriedyGroup.com)