

Accutron™ Nasal Masks Overview



HuFriedyGroup is the only manufacturer of three different designs of single-use nasal masks—so there is a style that will meet every clinician's needs!



Accutron™ ClearView™ Nasal Mask

The most popular Accutron design! This “mask-in-mask” style keeps the gas enclosed without the use of a scavenging hub and has an improved facial seal for added comfort and safety.



Accutron™ Axess™ Nasal Mask

The lowest profile Accutron design, with the lightest weight tubing. No cannula or tape needed, so it does not sacrifice comfort for fit!



Accutron™ PIP+™ Nasal Mask

The original Accutron™ design and still a reliable choice for nitrous delivery. Available with both single- and multi-use options.

	ClearView™	Axess™	PIP+™	PIP+™ Multi-Use
Single-Use	✓	✓	✓	
Individually Wrapped	✓	✓	✓	N/A
Visual Breath Monitoring	✓	✓		
Three Sizes	✓	✓	✓	✓
Fully enclosed/No Scavenging Hub	✓	✓		
Variety of Colors/Scents	✓	✓	✓	
Unscented Available	✓	✓	✓	✓
Uses Reservoir Bag/"Patient Demand" System	✓		✓	✓
Uses Continuous/Positive Flow System		✓		
Made in USA	✓	✓	✓	✓

Accutron is now a proud member of

HuFriedyGroup
The Best In Practice

Accutron™ Nasal Masks

Ordering Options



Accutron™ ClearView™ Classic Nasal Mask (12/box)	
33037-9	Pediatric Unscented Grey
33037-10	Pediatric Outlaw Orange
33037-11	Pediatric Sassy Strawberry
33037-12	TOP SELLER Pediatric Birthday Bubblegum
33037-17	Pediatric Groovy Grape
33035-9	Adult Unscented Grey
33035-10	Adult Outlaw Orange
33035-12	Adult Birthday Bubblegum
33035-14	TOP SELLER Adult French Vanilla
33035-16	Adult Fresh Mint
33034-9	Large Adult Unscented Grey
33034-10	Large Adult Outlaw Orange
33034-14	Large Adult French Vanilla
33034-16	Large Adult Fresh Mint

Accutron™ ClearView™ Nasal Masks connect effortlessly to:
 ClearView™ Scavenging Circuits I (**43005** and **43006**)
 ClearView™ Scavenging Circuit II (**43003** and **43004**)
 Any existing PIP+ Scavenging Circuit
 Any Porter™ Brown Scavenging Circuit (no adaptor required)
 Matrix Scavenging System (has one smooth tube on each side and comes with two connections) requires adaptor SKU **35360**

Accutron™ Axess™ Low Profile Nasal Mask (24/box)	
53037-9	Small Clear Unscented
53037-12	TOP SELLER Small Birthday Bubblegum
53037-16	Small Fresh Mint
53035-9	TOP SELLER Medium Clear Unscented
53035-12	Medium Birthday Bubblegum
53035-16	Medium Fresh Mint
53034-9	Large Clear Unscented
53034-12	Large Birthday Bubblegum
53034-16	Large Fresh Mint

Accutron™ PIP+™ Single-Use Nasal Masks (24/box)	
33017-1	TOP SELLER Small Variety Pack 1 (Orange, Strawberry and Bubblegum)
33017-9	Small Unscented Grey
33017-10	Small Outlaw Orange
33017-11	Small Sassy Strawberry
33017-12	Small Birthday Bubblegum
33017-17	Small Groovy Grape
33016-1	TOP SELLER Medium Variety Pack 1 (Orange, Strawberry and Bubblegum)
33016-2	Medium Variety Pack 2 (Pina Colada, Vanilla and Orange)
33016-3	Medium Variety Pack 3 (Mint, Peach and Vanilla)
33016-9	Medium Unscented Grey
33016-10	Medium Outlaw Orange
33016-11	Medium Sassy Strawberry
33016-12	Medium Birthday Bubblegum
33016-13	Medium Pina Colada
33016-14	Medium French Vanilla
33016-15	Medium Perfect Peach
33016-16	Medium Fresh Mint
33015-3	TOP SELLER Large Variety Pack 3 (Mint, Peach and Vanilla)
33015-9	Large Unscented Grey
33015-14	Large French Vanilla
33015-15	Large Perfect Peach
33015-16	Large Fresh Mint

Accutron™ PIP+™ Multi-Use Nasal Masks (Steam Sterilize) (1 each)	
33018	Multi-Use Large Mask
33019	Multi-Use Medium Mask
33020	Multi-Use Small Mask

For more information, visit [HuFriedyGroup.com/Accutron](https://www.HuFriedyGroup.com/Accutron)

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | [HuFriedyGroup.com](https://www.HuFriedyGroup.com)
 All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.
 ©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. 0622

Accutron is now a proud member of

