




HU-FRIEDY KEY OPINION LEADER

ACADEMY OF RICHARD V. TUCKER STUDY CLUBS KIT LIST

The Academy of Richard V Tucker Study clubs

consist of 75 independently operating, clinical study clubs throughout the United States, Canada, Germany, Italy, and Japan. Each study club consists of approximately 12 members and a mentor. The mentor ensures that the study club session adheres to the Tucker technique. Each session is divided into a clinical component and a critique.



PART CODE	DESCRIPTION	QUANTITY
MT232TM9	232 TruBal Margin Trimmer #9Hdl	
MT233TM9	233 TruBal Margin Trimmer #9Hdl	
MT232TBM9	232 Tucker Margin Trimmer #9Hdl	
MT233TBM9	233 Tucker Margin Trimmer #9Hdl	
CP14/14-09	#14/14-0 DE Hatchet (15-10-16) #9Hdl	
CP15/15-09	#15/15-0 DE Hatchet (15-10-16) #9Hdl	
CP44S9	#44S Off-Angle Cavity Prep #9Hdl	
CP45S9	#45S Off-Angle Cavity Prep #9Hdl	
IM5109	10 Inst. Sig Cassette Green	
KIT-TUCKER	Complete Tucker Kit	

Check One: Quote ☐ Order ☐

Contact 1-800-Hu-Friedy for more information or fax order to 1-800-729-1299.

Name _____

Address _____

City _____

State _____

Zip Code _____

Phone _____

Preferred Dealer _____

Dealer Rep Name _____

Please check the box if you would like to receive confirmation: ☐



Visit us online at Hu-Friedy.com

How the best perform



THE ACADEMY OF RICHARD V TUCKER STUDY CLUBS KIT LIST

A MT232TM9

B MT233TM9

C MT232TBM9

D MT233TBM9


E CP14/14-09

F CP15/15-09

G CP44S9

H CP45S9

I IM5109



THE ACADEMY OF RICHARD V TUCKER STUDY CLUBS KIT LIST



- A MT232TM9**
232 TruBal Margin Trimmer #9Hdl
Margin Trimmer used to produce proper bevel on enamel margins. Similar to a hatchet except the blade is curved and the cutting edge angled.



- B MT233TM9**
233 TruBal Margin Trimmer #9Hdl
Margin Trimmer used to produce proper bevel on enamel margins. Similar to a hatchet except the blade is curved and the cutting edge angled.



- C MT232TBM9**
232 Tucker Margin Trimmer #9Hdl
Margin Trimmer used to produce proper bevel on enamel margins. Similar to a hatchet except the blade is curved and the cutting edge angled. Margin Trimmer used to produce proper bevel on enamel margins. Similar to a hatchet except the blade is curved and the cutting edge angled.



- D MT233TBM9**
233 Tucker Margin Trimmer #9Hdl
Margin Trimmer used to produce proper bevel on enamel margins. Similar to a hatchet except the blade is curved and the cutting edge angled.



- E CP14/14-09**
#14/14-0 DE Hatchet (15-10-16) #9Hdl
Used for cavity preparation: retentive areas, internal line angles and removing hard caries. 1.5mm/1.5mm [15-10-16]



- F CP15/15-09**
#15/15-0 DE Hatchet (15-10-16) #9Hdl
Used for cavity preparation: retentive areas, internal line angles and removing hard caries. 1.5mm/1.5mm [15-10-16]



- G CP44S9**
#44S OFF-ANGLE CAVITY PREP #9HDL
Used for cavity preparation: retentive areas, internal line angles and removing hard caries.



- H CP45S9**
#45S Off-Angle Cavity Prep #9Hdl
Used for cavity preparation: retentive areas, internal line angles and removing hard caries.



- I IM5109**
10 Inst. Sig Cassette Green
IMS® cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.