

HuFriedyGroup Key Opinion Leader

Prof. Istvan Urban DMD, MD, PhD
General Surgery Kit & Sinus Lift Kit



Focus on: **Guided Bone
Regeneration**

Country: **Hungary**

Favourite Instrument:
Mini Me (PFIWDS1MKX)

Website:
implant.hu/en

Inspired by: **Henry H. Takei**

VISIT THE
KOL COMMUNITY



I enjoy working without a compromise.

Prof. Istvan Urban DMD, MD, PhD




GENERAL SURGERY KIT

Instruments are placed in a surgical cassette according to the order in which they will be used during the surgery. These instruments are sterilized in the surgical cassette before the surgery. The instruments are used out of the cassette and are kept in order during the entire surgical procedure.


URBANGSKIT

Istvan Urban General Surgery Kit



GENERAL SURGERY KIT

Instruments Designed by Istvan Urban



The Urban Mini Me

is designed with angled, sharp working ends for minimal invasive flap elevation in the esthetic zone and for papilla elevation and tunneling. The Mini Me can be used for much more—for lingual flap elevation, used during mandibular ridge augmentation surgery and for blunt flap advancement.

9



PFIWDS1MKX
Mini Me Microsurgical
Periosteal, **Black Line**

7



CURBAN1
Urban Chisel

6



SC4R/4L9E2
4R/4L EverEdge™ 2.0 Columbia Curette,
#9 Handle

8



SB2/36
B-2/3 Ball DE Scaler, #6 Handle

2



MIR5HD
#5 HD Single-Sided Mirror

1

MH6
#6 Mirror Handle, Cone Socket

GENERAL SURGERY KIT



3
PCNC15/126
UNC 15/12 Color Probe,
#6 Handle



4
10-130-05E
5 Scalpel Handle



5
PT1X
1 Periotome, Posterior,
Black Line




10
PPR3X
PR-3 Prichard Periosteal,
Short, **Black Line**

Black Line

is designed to optimize clinical outcomes by delivering efficiency throughout perio and surgical procedures. Black Line instruments feature a dark performance engineered coating which reduces light reflection. The dark coating provides enhanced contrast, improving visual acuity at the surgical site, supra- and sub- gingivally. The coating also assists with edge retention.



13
CRM
University of Minnesota
Cheek Retractor




SuperCut
Cutting blade honed/
serrated to "razor"
sharpness with opposing
blade finished with fine
serrations to avoid
tissue slippage.

17
S16SC
16 Goldman-Fox Scissors, SuperCut,
Curved, Serrated, 12.5 cm

GENERAL SURGERY KIT



11



TPSLCOSM


Micro Cooley Tissue Pliers,
Straight, Modified, 18 cm

14



NH5024SIM

5024R Simion **Perma Sharp™** Needle Holder,
Round, Straight, 18 cm



12



TP5090

DeBakey **Perma Sharp™** Tissue Pliers,
Straight, 15 cm

Perma Sharp™ Instruments

Featuring tungsten carbide inserts bonded to surgical-grade stainless steel on both cutting edges for smooth, controlled cutting and maximum durability.

Did you know?

All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like Instrument Lubricant Spray (**ILS**), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilized in the open position.**



18



H4

4 Halsted-Mosquito Hemostat,
Straight, 12 cm

15



DP17

17 Dressing Pliers

GENERAL SURGERY KIT

16

IMS-1372SH

IMS™ 1/2 Bur
Cushion, Short
Lid, 6 Burs



19

MP54

Markel Mouth Prop,
Adult



Additional instruments not included in the General Surgery Kit



ASPCG3

Suction Tube Cogswell



CM2

2 Molt Surgical Curette, #502 Handle

WHY YOU SHOULD USE AN INSTRUMENT MANAGEMENT SYSTEM:

By standardizing your instruments with IMS, additional time and revenue is generated while extending the useful life of your HuFriedyGroup instruments. Dental professionals state that by implementing IMS™ Cassettes, they save 5 to 10 minutes per procedure allowing the staff to spend more time on value-added activities*. In addition, safer transportation of soiled instruments reduces risk of injury and exposure to blood borne pathogens.

* When compared to single instrument decontamination.
Based on market survey results.

DISCOVER
MORE AT **IMS™**
WEBPAGE




20



IM5101

16-Instrument IMS™ Infinity Series™
Cassette, Green, 4 Hinge Clips

SINUS LIFT KIT



URBANSLKIT

Istvan Urban Sinus Lift Kit

Instruments Designed by Istvan Urban



IMPURBAN16

Urban Sinus Lift Instrument 1, #6 Handle



IMPURBAN26

Urban Sinus Lift Instrument 2, #6 Handle



IMPURBAN36

Urban Sinus Lift Instrument 3, #6 Handle

SINUS LIFT KIT

5



PPR3S6

3 Prichard Periosteal, Short #6 Handle

1



PLGJO/46

4 Jovanovic Plugger, #6 Handle

6



CRURBAN

Urban Cheek Retractor

7



IM5101

8-Instrument IMS™ **Infinity Series™** Cassette
w/ Accessory Area Dimensions 152 x 184 x 34 mm

Additional instruments not included in the Sinus Lift Kit

ASPON

Pontius Suction Tube



URBAN REGENERATION INSTITUTE

COURSES:

- 9-Day Regenerator Program
- 3-Day Advanced Bone and Soft Tissue Course
- 2-Day Soft Tissue Course
- 2-Day Sinus Course
- 4-Day IDR Course
- Endo Course



Developing minimally invasive regenerative techniques led to favorable long-term clinical outcomes. I have a passion for sharing those important surgical details, the science behind all the innovation and the long term follow up of the cases that puts all details to a credible prospective.

I would like to invite you to join our courses where we honestly discuss the realities of hard and soft tissue regeneration.

Prof. Istvan Urban



WATCH
THE COURSES
VIDEO



For more information visit [HuFriedyGroup.com](https://www.HuFriedyGroup.com)