

TECHNICAL DATA SHEET

PRODUCT INFORMATION

Sparkle™ Prophy Paste with Xylitol

SKU	Description	Quantity
UPC12	Course grit, Mint	12 oz/340 g jar
UPCB	Course grit, Berrylicious Assortment*	200/Box
UPCBG	Course grit, Bubble Gum	200/Box
UPCCH	Course grit, Cherry	200/Box
UPCM	Course grit, Mint	200/Box
UPCOV	Course grit, Orange Vanilla	200/Box
UPCA	Course grit, Assorted Flavors	200/Box
UPMB	Medium grit, Berrylicious Assortment*	200/Box
UPMBG	Medium grit, Bubble Gum	200/Box
UPMCH	Medium grit, Cherry	200/Box
UPMM	Medium grit, Mint	200/Box
UPMOV	Medium grit, Orange Vanilla	200/Box
UPMA	Medium grit, Assorted Flavors	200/Box
UPFBG	Fine grit, Bubble Gum	200/Box



PRODUCT OVERVIEW

Sparkle™ Prophy Paste is a blend of polishing and cleaning agents, which makes it effective in removing stains and smoothing the surfaces of teeth. Variety of available grits allows you to choose the least abrasive polishing agent to remove biofilm and stain.

MATERIALS

- Glycerin
- Water
- Pumice
- Sodium Fluoride
- FD&C Coloring and flavoring
- Xylitol

DETAILS

- Contains 1.23% Fluoride ion
- Formulated to maximize stain removal with minimal enamel loss
- Advanced formulation provides consistent, splatter-free results
- Completely rinses out, leaves no residual grit
- Available in multiple grits and flavors

TECHNICAL DATA SHEET

ALLERGEN INFORMATION

Product is free of gluten, soybeans/soy products, milk/milk products, peanuts/peanut products, tree nuts/tree nut derivatives, hazelnuts, egg/egg products, banana, pineapple, wheat/wheat products, fish/fish products, shellfish, and sesame

STORAGE REQUIREMENTS

- n/a

SHELF LIFE

3 years

REGULATORY INFORMATION

Regulation: 872.6030, Oral Cavity Abrasive

Polishing Agent

Medical Device Listing: D445899

Product Code: EJR Classification: FDA Class 1

*Includes: Raspberry, Blueberry, Strawberry

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | HuFriedyGroup.com

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.

©2022 Hu-Friedy Mfg. Co., LLC. All rights reserved.