

TECHNICAL DATA SHEET

PRODUCT INFORMATION

Sparkle™ FREE Prophy Paste with Xylitol

SKU	Description	Quantity
UPSFCC	Course grit, Cinnamon	200/Box
UPSFCE	Course grit, Fruity	200/Box
UPSFCF	Course grit, Spearmint	200/Box
UPSFCWC	Course grit, White Chocolate	200/Box
UPSFMC	Medium grit, Cinnamon	200/Box
UPSFME	Medium grit, Fruity	200/Box
UPSFMS	Medium grit, Spearmint	200/Box
UPSFMWC	Medium grit, White Chocolate	200/Box



PRODUCT OVERVIEW

Sparkle™ Free Prophylaxis Paste is a blend of polishing and cleaning agents, which makes it effective in removing stains and smoothing the surfaces of teeth. Variety of available grits allows you to choose the least abrasive polishing agent to remove biofilm and stain.

MATERIALS

- Glycerin
- Pumice
- Water
- Xylitol

DETAILS

- Prophy Paste in a unique fluoride/dye free formula
- Formulated to maximize stain removal with minimal enamel loss
- Advanced formulation provides consistent, splatter-free results
- Completely rinses out, leaves no residual grit
- Available in multiple grits and flavors

ALLERGEN INFORMATION

- Product is free of nuts, soy, milk, eggs, gluten, aspartame and corn/corn products and is manufactured in area in which has no exposure to the following items: wheat, eggs, peanut, milk, fish, nut, soybean and shellfish

STORAGE REQUIREMENTS

- n/a

TECHNICAL DATA SHEET

SHELF LIFE

3 years

REGULATORY INFORMATION

Product Code: EJR

Classification: FDA Class 1

Regulation: 872.6030, Oral Cavity

Abrasive Polishing Agent

Medical Device Listing: D445899