

# TECHNICAL DATA SHEET

## PRODUCT INFORMATION

### Sparkle™ FREE Prophy Paste with Xylitol

SKU	Description	Quantity
UPSFCC	Course grit, Cinnamon	200/Box
UPSFCE	Course grit, Fruity	200/Box
UPSFCS	Course grit, Spearmint	200/Box
UPSFCEW	Course grit, White Chocolate	200/Box
UPSFMC	Medium grit, Cinnamon	200/Box
UPSFME	Medium grit, Fruity	200/Box
UPSFMS	Medium grit, Spearmint	200/Box
UPSFMEW	Medium grit, White Chocolate	200/Box



## PRODUCT OVERVIEW

Sparkle™ Free Prophylaxis Paste is a blend of polishing and cleaning agents, which makes it effective in removing stains and smoothing the surfaces of teeth. Variety of available grits allows you to choose the least abrasive polishing agent to remove biofilm and stain.

## MATERIALS

- Glycerin
- Pumice
- Water
- Xylitol

## DETAILS

- Prophy Paste in a unique fluoride/dye free formula
- Formulated to maximize stain removal with minimal enamel loss
- Advanced formulation provides consistent, splatter-free results
- Completely rinses out, leaves no residual grit
- Available in multiple grits and flavors

## ALLERGEN INFORMATION

- Product is free of gluten, soybeans/soy products, milk/milk products, peanuts/peanut products, tree nuts/tree nut derivatives, hazelnuts, egg/egg products, banana, pineapple, wheat/wheat products, fish/fish products, shellfish, and sesame

## STORAGE REQUIREMENTS

- n/a

# TECHNICAL DATA SHEET

## SHELF LIFE

---

3 years

## REGULATORY INFORMATION

---

Product Code: EJR

Classification: FDA Class 1

Regulation: 872.6030, Oral Cavity

Abrasive Polishing Agent

Medical Device Listing: D445899