


VistaPure™

Water Purification System



Meet your steri-center's new best friend.

The **VistaPure™ Water Purification System** provides two grades of high-quality water for steam sterilizers, dental bottles, ultrasonic cleaners and the final-rinse cycle of compatible instrument washers. The system connects directly to the city water supply and sends high-quality water anywhere you need it in the steri-center for easy, on-demand use.



Two grades of the water practices need

For steam sterilizers:

Point-and-shoot wand makes filling sterilizers with autoclave-ready water a breeze. No more lifting heavy jugs of water!

For bottle systems:

Mounted faucet provides non-corrosive bottle-quality (hyperfiltered) water for bottle systems. No more buying bottled water!



Instant Status Monitoring

The VistaPure™ Water Purification System features an inline TDS meter that monitors water quality — and makes it easy to know when it's time to replace your filters.



No Electricity Required

The VistaPure™ Water Purification System operates on water and air pressure, so no electricity is required to supply all the high-purity water your office needs!



Environmentally Friendly

The VistaPure™ Water Purification System doesn't require chemicals, and its efficient design minimizes the amount of wastewater produced during filtration.

Crosstex is now a proud member of

HuFriedyGroup
The Best in Practice

VistaPure™

Water Purification System



Optional Accessories

Instrument Washer Adapter Kit CVS3010

Protect expensive instruments and equipment by using the VistaPure™ Water Purification System's hyperfiltered water for the final-rinse cycle of instrument washers like the HYDRIM® Automated Instrument Washer.

AutoFill for BRAVO™ Autoclave Chamber CVS3030


Replenish the SciCan BRAVO™ Autoclave Chamber's refill tank with autoclave-quality water automatically — enabling truly hands-free convenience!

AutoFill for STATIM® Autoclave CVS3020

Deliver autoclave-quality water to your SciCan STATIM® Autoclave automatically. You'll never have to fill your STATIM Autoclave manually again!



Space Requirements



5-Year Filter-Replacement Projection *

INSTALLATION	YEAR 2	YEAR 3	YEAR 4	YEAR 5	TOTALS (USD)
1	1	1	1	1	\$321
2		2	2	2	\$348
3		3	3	3	\$348
4				4	\$232
5	5	5	5	5	\$321
					\$1570
					Average of \$314 per year

Learn More at [HuFriedyGroup.com](https://www.hufriedygroup.com)

*Assumes deionization filters are replaced every 18 months and hyperfiltration element is replaced every four years, on average. Ultimately, actual filter lifespan depends upon local water quality, inbound water pressure and volume of water used.

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | [HuFriedyGroup.com](https://www.hufriedygroup.com)

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC its affiliates or related companies, unless otherwise noted. BRAVO™, HYDRIM® and STATIM® are trademarks and registered trademarks of SciCan Corporation Ltd.

©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. DLIT00637 Rev D 1022

Crosstex is now a proud member of
HuFriedyGroup
The Best in Practice