


Diamond Dusted Forceps

Our **Diamond Dusted Forceps** feature a fine diamond coating which increases resistance between the forceps and tooth surface providing enhanced engagement. Diamond Dusted Forceps are available in an array of patterns which allows for full-mouth coverage.

A light diamond coating bonded to longitudinal beak serrations increases grip by 20% compared to traditional forceps.*



Thin, tapered beak designs allow for improved access compared to traditional forceps beaks.



Matte finish for additional contrast and reduced light reflection.

Unique hole design reduces the weight of the forceps, providing increased comfort and control.

Diamond Dusted Forceps

are available in the following Atraumair patterns:



Upper

F1XSD **1 Standard Forceps, Diamond Dusted**
Upper Incisors



FAF1XSD **1 Apical Forceps, Diamond Dusted**
Upper Incisors




F35XSD **35 Forceps, Diamond Dusted**
Upper Canines and Premolars




FAF35XSD **35 Apical Forceps, Diamond Dusted**
Upper Canines and Premolars




FX44XSD **44 Root Forceps, Diamond Dusted**
Upper Roots



FX51XSD **51 Root Forceps, Diamond Dusted**
Upper Roots



F67XSD **67 Apical Forceps, Diamond Dusted**
Upper Molars



Lower

F13XSD **13 Forceps, Diamond Dusted**
Lower Premolars



F21XSD **21 Forceps, Diamond Dusted**
1st and 2nd Molars




F36XSD **36 Forceps, Diamond Dusted**
Lower Incisors



FAF36XSD **36 Apical Forceps, Diamond Dusted**
Lower Premolars and Incisors




FX45XSD **45 Root Forceps, Diamond Dusted**
Lower Roots, Serrated



FX73XSD **73 Forceps, Diamond Dusted**
Lower Molars



FX79XSD **79 Forceps, Diamond Dusted**
Lower Molars



Find your optimal extraction forceps at [HuFriedyGroup.com/Surgical-Forceps](https://www.hufriedygroup.com/Surgical-Forceps)

*Data on file compared to standard Hu-Friedy Forceps.

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | [HuFriedyGroup.com](https://www.hufriedygroup.com)

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC its affiliates or related companies, unless otherwise noted.

©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-666/0123

