

Key Opinion Leader

Dr. Paulo Maló

ALL-ON-4 & ZYGOMA KIT

A full solution for a variety of implant surgeries



VISIT THE
KOL COMMUNITY



Focus on: **Oral and
Maxillofacial Surgery**



Country: **Portugal**



Favorite Instrument:
Malo Surgical Retractor



Website:
www.malodental.com



Inspired by: **Prof. Brånemark**



This surgical setup does not represent only a complete solution to perform the "All-on-4" technique but in fact it can also be used for all kinds of implants surgeries and of course zygomatic procedures - or as I call it "All-on-4 hybrid" - using the two additional and specifically designed instruments.

Dr. Paulo Maló



TWO NEW INSTRUMENTS FOR ZYGOMATIC PROCEDURES

The instrumentation for "All-on-4" procedures, identified by Dr. Maló, has been extended and expands to an even wider field - the zygomatic approach. The surgical retractor (SRMALO) is sharp and curved and helps to push away the soft tissue of the zygoma. The long handled periosteal (PMALO) provides better reach, improving field of vision. Both instruments are a crucial addition to the "All-on-4" surgical kit, allowing for zygomatic procedures.




SRMALO

Malo Surgical Retractor



PMALO

Malo Periosteal



MALALL4KIT

All-on-4 Kit



WATCH
THE **NEW**
INSTRUMENTS
VIDEO

WATCH
THE **ALL ON 4**
VIDEO



Visit us online at **HuFriedyGroup.com**

HuFriedyGroup
The Best In Practice

A FULL SOLUTION FOR ALL KINDS OF IMPLANT SURGERIES

CHOOSE THE COLOR.

Just add the color coding number to the last part of the partcode (replace “_” with the number)

Red	●	0
Grey	●	1
White	●	2
Orange	●	3
Purple	●	4
Yellow	●	5
Lavender	●	6
Black	●	7
Blue	●	8
Green	●	9
Ocean Blue	●	11

WHY YOU SHOULD USE AN INSTRUMENT MANAGEMENT SYSTEM (IMS):

By standardizing your instruments with IMS™ System, additional time and revenue is generated while extending the useful life of your HuFriedyGroup instruments. Dental professionals state that by implementing IMS™ Cassettes, they save 5 to 10 minutes per procedure allowing the staff to spend more time on value-added activities*. In addition, safer transportation of soiled instruments reduces risk of injury and exposure to blood borne pathogens.


*When compared to single instrument decontamination. Based on market survey results.

DISCOVER
MORE AT IMS™
WEBPAGE



MALALL4KIT

All-on-4 Kit (IMS™ Infinity Series
Cassette red, included)



LEARN MORE
ABOUT EACH
INSTRUMENT BY
CLICKING ON THE
ICON NEXT TO THE
DESCRIPTION



NHLBPSB

Needle Holder Lichtenberg
PermaSharp™, 20cm



TP5091

Tissue Pliers De Bakey straight
PermaSharp™, 19cm



10-130-05

#5 Scalpel Handle, Round Straight



CRB

Cheek Retractor Bishop



CL846

Surgical Curette, Lucas #84
2.2mm



CL886

Surgical Curette, Lucas #86
2.8mm



PP5689

Malo Periosteal



DID YOU KNOW?

All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilised in the open position.**



XP3A/NC156

Expro, North Carolina #3A
15-14-13....3-2-1



S5002

Scissors, Kelly Straight
PermaSharp™ 16cm



RB EY

Rongeur Beyer, 18cm



IR

Retractor Vestibule Hook
Universal, 20cm





PMALO



SRMALO

DO YOU WANT TO **EXPAND** THE ALL-ON-4 KIT INTO A FULL ZYGOMA SETUP? MAKE SURE YOU **ADD** THE NEW **MALO PERIOSTEAL** AND **MALO SURGICAL RETRACTOR**!

WATCH
THE **NEW**
INSTRUMENTS
VIDEO



LEARN MORE ABOUT **MALO DENTAL A GLOBAL DENTAL COMPANY**



Focus on: **Oral and Maxillofacial Surgery**



Country: **Portugal**



Favorite Instrument:
Malo Surgical Retractor



Website:
www.malodental.com



Inspired by: **Prof. Brånemark**



MALO DENTAL is born, bringing together many colleagues globally, under my leadership. My aim has been to ensure that the MALO brand would continue, as it is the only dental brand which has international presence, being represented in 16 countries at the moment, encompassing more than 150 clinics either directly established as Malo Dental, or in a partnership with Malo Dental. The MALO brand is synonymous with globally recognised credibility, numerous publications, and has developed products and techniques utilised internationally. The MALO brand is without question a highly successful international brand, which has already been used as a case study. I have dedicated myself for more than 25 years to developing techniques and products which have contributed to the recognition I have, and have learnt itself to the prestige of the MALO brand.

Dr. Paulo Maló



For more information visit **HuFriedyGroup.com**