

Surgical Kit

RESCUE Kit (IMRESCUE)



1 IM4160
16 Instrument Signature Series™
Cassette, Red

IMS™ Cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.



2 IMS-1372SH
IMS 1/2 Bur Cushion - Short Lid

Designed to hold burs securely in place during cleaning and sterilization. Hollow design allows drainage and proper cleaning.



3 EPTSSP
Luxating Hybrid Elevator, Straight Spade

A unique cross between a periosteal and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



4 EPTSMC
Luxating Hybrid, Small Curved

A unique cross between a periosteal and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



5 EPTSMS
Luxating Hybrid Elevator, Small Straight

A unique cross between a periosteal and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



6 CL866
86 Lucas DE Surgical Curette #6 Handle

Designed for curettage, cyst removal and tooth socket debridement. 2.8mm.



7 EHB13/146
13/14 Heidbrink DE Root Pick


Used for the removal of root tips in difficult-to-reach areas. Designed for delicate root teasing procedures.



8 RMF
Mini Friedman Rongeurs

Used for trimming and recontouring alveolar bone and gross tissue removal.

Rescue Kit (IMRESCUE)



Not pictured: #2 [IMS-1372SH](#)

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | [HuFriedyGroup.com](#)
All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.
©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-351/1023