

Surgical Kit

Rescue Kit, Black Line (IMNRESCUEX)



- 1** **IMN8120**
12 Instrument Infinity Series™ Cassette, Red
IMS™ cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.



- 2** **EPTSMSX**
Luxating Hybrid Elevator, Small Straight, Black Line
A unique cross between a periotome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



- 3** **EPTSMCX**
Luxating Hybrid, Small Curved, Black Line
A unique cross between a periotome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



- 4** **EPTSSPX**
Luxating Hybrid Elevator, Straight Spade, Black Line
A unique cross between a periotome and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.



- 5** **EHB13/14X**
13/14 Heidbrink DE Root Pick, Black Line
Used for the removal of root tips in difficult-to-reach areas.




- 6** **CL86X**
86 Lucas DE Surgical Curette #6 Handle, Black Line
Designed for curettage, cyst removal and tooth socket debridement. 2.8mm



- 7** **RMF**
Mini Friedman Rongeurs
Used for trimming and recontouring alveolar bone and gross tissue removal. 1.5mm/14.5cm

Rescue Kit, Black Line (IMNRESCUEX)



Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | HuFriedyGroup.com
All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.
©2023 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-333/1023