

Surgical Kit

Lateral Sinus Lift Kit, Black Line (IMNLATKITX)



- 1** **IM4162**
16 Instrument Infinity Series™ Cassette, White
IMS™ Cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.



- 2** **IMP6578SX**
Kramer-Nevins Sinus Lift Instrument, 47-Degree, Black Line
Used during a sinus lift procedure to separate/reflect the Schneiderian membrane from the maxillary bone and to elevate the membrane. With obtuse angles 4.1mm/7.6mm.



- 3** **IMP6577SCX**
Kramer-Nevins Sinus Lift Instrument, Strong Curve, Black Line
Used during a sinus lift procedure to separate /reflect the Schneiderian membrane from the maxillary bone and to elevate the membrane. With acute angles 4.1mm/7.6mm.



- 4** **BPACK**
Bone Packer DE 6-8mm
Used to condense bone during grafting procedures with working ends 6mm and 8mm wide.



- 5** **SYBG**
Syringe Bone Grafting 7mm
Used during grafting procedures. Inner Diameter 5.8 mm. Load from tip to prevent clogging.



- 6** **BONT**
Bone Tamper 7mm
Used to separate the Schneiderian membrane from surrounding bone during lateral sinus lift procedures.




- 7** **NHCW**
Crile-Wood Needle Holder
Used to grasp and guide the needle when suturing. Recommended suture sizes: 3-0, 4-0, 5-0, 6-0



- 8** **S16SCX**
16 Goldman-Fox Scissors, Black Line
Used for trimming tissue or cutting sutures. Featuring a non-reflective surface, improved contrast between tissue/scissor. 13.5cm.




- 9** **CRURBAN**
Urban Cheek Retractor
Used to hold mucoperiosteal flaps, cheeks, lips and tongue away from the surgical area. 20.6mm



- 10** **IM1002**
6-Pack Hinged Instrument Clip
Secure and maintain open position of scissors, hemostats, and other hinged instruments in lid of cassette.

Lateral Sinus Lift Kit, Black Line (IMNLATKITX)



Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | HuFriedyGroup.com
All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.
©2023 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-348 /1023