

Surgical Kit

Crestal Sinus Lift Kit, Black Line (IMNCREKITX)



1 IMN4162
16 Instrument Infinity Series™
Cassette, White

IMS™ Cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.



2 10-130-03
3 Scalpel Handle

Used for cutting gingival tissue and making surgical incisions. Handle is 4.9" long and 13mm wide at base. Compatible with standard scalpel blade sizes.



3 P9X
9 Molt Periosteal Elevator, Black Line

Designed for reflecting and retracting the mucoperiosteum after incisions of the gingival tissue. 4.3mm/8.4mm



4 CM2/4X
2/4 Molt DE Curette, Black Line

Designed for curettage, cyst removal and tooth socket debridement with 3.6mm and 7mm working ends.



5 TP5043
Adson-Brown Perma Sharp™
Tissue Pliers

Used to grasp and stabilize soft tissue during suturing. 12cm



6 CSEUWES6
Euwe Bone Scoop/Plugger, #6 Handle

Designed for easy transport and non-stick compaction of bone graft material. 3.0mm/7.0mm.



7 PCP65006
20mm Probe 8-10-12-14-16-18-20,
#6 Handle

A double-ended instrument with an explorer on one end and a probe on the other, excellent for dental examinations and charting.



8 OSTPAL1
1 Palti Osteotome 2.7mm

Used for preparing bone. 7-10-13-15-18mm



9 OSTPAL2
2 Palti Osteotome 3.2mm

Used for preparing bone. 7-10-13-15-18mm



10 OSTPAL3
3 Palti Osteotome 3.7mm

Used for preparing bone. 7-10-13-15-18mm



11 OSTPAL4
4 Palti Osteotome 4.2mm

Used for preparing bone. 7-10-13-15-18mm



12 PFIWDS16
1 Woodson Plastic Filling, #6 Handle


Designed to reflect and retract the mucoperiosteum after incisions of the gingival tissue with 2.2mm and 3.8mm width.



13 NH5038
Crile-Wood Perma Sharp™ Needle Holder

Used to grasp and guide the needle when suturing. Featuring 3-0 to 5-0 carbide beaks, finger ring style.

Crestal Sinus Lift Kit, Black Line (IMNCREKITX)



Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | HuFriedyGroup.com
All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.
©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-3458/1023