



HU-FRIEDY PROCEDURAL SET-UP

BLACK LINE RUBBER DAM CLAMP KIT



PART CODE	DESCRIPTION	QUANTITY
RDCM13AX	Rubber Dam Clamp #13A, Black Line	
RDCM12AX	Rubber Dam Clamp #12A, Black Line	
RDCM14X	Rubber Dam Clamp #14, Black Line	
RDCM27NX	Rubber Dam Clamp #27N, Black Line	
RDCMW8AX	Rubber Dam Clamp #W8A, Black Line	
RDCM212X	Rubber Dam Clamp #212, Black Line	
RDCM2AX	Rubber Dam Clamp #2A, Black Line	
RDCM26X	Rubber Dam Clamp #26, Black Line	
RDCOB	Rubber Dam Clamp Organizing Board	
RDCKITX	Complete Black Line Rubber Dam Clamp Kit	

BLACK LINE RUBBER DAM CLAMP KIT



- A RDCM13AX**
Rubber Dam Clamp #13A, Black Line
Serrated offset jaws for lower left molars; paired with RDCM12AX.



- B RDCM12AX**
Rubber Dam Clamp #12A, Black Line
Serrated offset jaws for lower right molars; paired with RDCM13X.



- C RDCM14X**
Rubber Dam Clamp #14, Black Line
Deeply festooned jaws for partially erupted or irregularly shaped molars.



- D RDCM27NX**
Rubber Dam Clamp #27N, Black Line
Slightly festooned jaws for premolars.



- E RDCMW8AX**
Rubber Dam Clamp #W8A, Black Line
Slightly festooned jaws for molars.



- F RDCM212X**
Rubber Dam Clamp #212, Black Line
Flat jaws for labial caries on anterior teeth.



- G RDCM2AX**
Rubber Dam Clamp #2A, Black Line
Flat jaws for larger premolars.



- H RDCM26X**
Rubber Dam Clamp #26, Black Line
Slightly reversed festooned jaws for small lower molars.



- I RDCOB**
Rubber Dam Clamp Organizing Board
Rubber Dam Clamp Organizing Board holds up to 8 clamps during sterilization and fits into cassette accessory area.