

HuFriedyGroup

Key Opinion Leader

Dr. Sascha A. Jovanovic

The Implant Surgery & Tissue Regeneration Kit



Focus on:
**Periodontal & Implant
Surgery**



Country: **US - California**



Favourite Instrument: **PJO/16
JO 1 - Jovanovic 4 mm Straight/
Backaction Chisel**



Website: **www.gidedental.com**



Inspired by:
Prof. Hubertus Spiekermann



University:
**UCLA, LOMA LINDA
& ACTA**



VISIT THE
KOL COMMUNITY



SJOVANKIT

Full Jovanovic Implant
Surgery & Tissue
Regeneration Kit




HuFriedyGroup instruments are the best high quality, precise hand instruments in the industry, and a natural extension of my hands.

Dr. Sascha A. Jovanovic



THE IMPLANT SURGERY & TISSUE REGENERATION KIT

The Jovanovic Implant & Regeneration Surgical Kit includes exceptionally crafted instruments specifically designed for advanced & minimal invasive periodontal and implant surgery. This sophisticated kit includes precision-engineered instruments fashioned from premium materials. Ergonomics plays a critical role in a clinician's long-term success which led to the incorporation of Black Line instruments. The dark coating minimizes glare, is visible when used subgingivally, and adheres less to surgical materials. The use of IMS™ Cassettes ensures easy access to instruments, keeping everything needed for procedures within reach.



Visit us online at **HuFriedyGroup.com**

THE IMPLANT SURGERY & TISSUE REGENERATION KIT



18

10-130-05E
Scalpel Handle
#5E European Style
Round, Matte Finish



19

SHDPV
Scalpel Handle Steel
Microsurgical



12

NHDPV
Diamond Dusted
Castroviejo
Microsurgical Needle
Holder, Straight, 18 cm



13

TPSLCOSM
Micro Tissue Plier,
Straight, 17.78 cm



11

PFIWDS1MKX
Microsurgical Periosteal
Black Line



17

MIR5HDB (2x)
#5 HD **Black Line**
Single-Sided Mirror



9

CRMX
University of Minnesota
Retractor, **Black Line**

16

MH6B (2x)
#6 **Black Line** Cone
Socket Mirror Handle

BLACK LINE
Black Line is designed to optimize clinical outcomes by delivering efficiency throughout the entire perio and surgical procedure. Black Line features a performance engineered coating for optimal edge retention, reduced light reflection afforded by the black matte finish and enhanced contrast and visual acuity at the surgical site and underlying tissue.



22

PESIM16
Sinus Lift Instrument
SE #6



4

CL86X
#86 Lucas Spoon
Shape Surgical
Curette, DE
Black Line



3

SKRA36
#KRA-3 Kramer
Periodontal Surgical
Curette DE,
#6 Handle



23

SGCSUB-06
Sub-0 Vision
Curvette #6



24

SC4R/4L9E2
Columbia 4R/4L
Curette #9



10

CRMS
Simion Cheek
Retractor



14

S14SCX
14 LaGrange Scissors,
Curved, Super-Cut
11.5 cm (4.5"), 1 serrated blade
Black Line




28

S20SCX
Goldman-Fox Scissors, Straight, Super-Cut
13.5 cm (5.25"), 1 serrated blade
Black Line


SUPER-CUT
Cutting blade honed/
serrated to "razor"
sharpness with opposing
blade finished with fine
serrations to avoid tissue
slippage.

THE IMPLANT SURGERY & TISSUE REGENERATION KIT


INSTRUMENTS DESIGNED IN COLLABORATION WITH DR. SASCHA A. JOVANOVIC



- PJO/16**
JO 1 - Jovanovic
4 mm Straight/
Backaction Chisel
#6 handle
- PJO/2M6**
JO 2 - Mini Jovanovic
2 mm Straight/
Backaction Chisel
#6 handle
- PJO/3M6**
JO 3 - Mini Jovanovic
3 mm Straight/
Backaction Chisel
#6 handle
- PLGJO/4X**
Jovanovic Bone Graft
Condenser DE
3 mm, Black Line
- PLGJO/56**
Jovanovic Bone Graft
Carrier/Condenser
#6 handle



- IMPJO66**
Jovanovic Sinus
Lift Instrument
#6 Handle
- IMPJO76**
Jovanovic Sinus
Lift Instrument #7,
#6 handle



- TP5090**
DeBakey Perma Sharp™
Tissue Plier with teeth,
15.24 cm
- TP5070**
General Perma Sharp™
Tissue Plier no teeth,
15.24 cm
- NH5024R**
Perma Sharp™ Castroviejo
Needle Holder, Round, 18 cm
- PCPUNC15CH**
Qulix SE UNC Probe, 1-15 Markings, #Smooth Handle
- PPR3X**
PR-3 Prichard
Periosteal, **Black Line**

**PERMA SHARP™
INSTRUMENTS**

Featuring tungsten carbide inserts bonded to surgical-grade stainless steel on both cutting edges for smooth, controlled cutting and maximum durability.

**CREATE YOUR
CUSTOM KIT AND
CHOOSE THE
CASSETTE WITH
YOUR FAVORITE
COLOR.**


Just add the color coding number to the last part of the partcode (replace " _ " with the number)

Red	0
Grey	1
White	2
Orange	3
Purple	4
Yellow	5
Lavender	6
Black	7
Blue	8
Green	9
Ocean Blue	11

WHY YOU SHOULD USE AN INSTRUMENT MANAGEMENT SYSTEM (IMS):

By standardising your instruments with IMS System, additional time and revenue is generated while extending the useful life of your HuFriedyGroup instruments. Dental professionals state that by implementing IMS™ Cassettes, they save 5 to 10 minutes per procedure allowing the staff to spend more time on value-added activities*. In addition, safer transportation of soiled instruments reduces risk of injury and exposure to blood borne pathogens.

*When compared to single instrument decontamination. Based on market survey results.



- IMN4203**
IMS™ Infinity Series™ **20 Instrument Cassette, Orange**
Dimensions: 203 mm x 280 mm x 32 mm

Also available and designed by Dr. Jovanovic:

JOVANOVIC ESSENTIAL IMPLANT KIT | SJOVANEKIT

This kit is designed to provide all of the essential instruments for Implant/Bone/Soft Tissue procedures.



PCPUNC15CH
Quilix SE UNC Probe,
1-15 Markings,
#Smooth Handle

10-130-05E
Scalpel Handle
#5E European Style
Round, Matte Finish

SHDPV
Scalpel Handle Steel
Microsurgical

PPR3X
PR-3 Prichard
Periosteal,
Black Line


PFIWDS1MKX
Microsurgical
Periosteal
Black Line



PESIM16
Sinus Lift Instrument SE #6

IMPJO76
Jovanovic Sinus Lift
Instrument #7, #6 handle

CL86X
#86 Lucas Spoon Shape
Surgical Curette, DE
Black Line



SKRA36
#KRA-3 Kramer Periodontal
Surgical Curette DE, #6 Handle

SC4R/4L9E2
Columbia 4R/4L Curette #9

PJO/16
JO 1- Jovanovic 4 mm Straight/
Backaction Chisel, #6 handle


PLGJO/56
Jovanovic Bone Condenser/
Periosteal, #6 handle

TP5090
DeBakey Perma Sharp™ Tissue Plier with teeth,
15.24 cm

NH5024R
Perma Sharp™ Castroviejo
Needle Holder, Round, 18 cm



S20SCX
Goldman-Fox Scissors, Straight, Super-Cut
13.5 cm (5.25"), 1 serrated blade
Black Line



IM1000
Hinged Instrument Clips
2-pack



IMN4163
Large Cassette With Accessory Area, Orange
Dimansions: 8" X 11" X 1.25" (203 mm x 280 mm x 32 mm)

SUGGESTED INSTRUMENTS (not included in the kit)



CM2/4X
2/4 Molt DE Curette **Black Line**



SG11/12R9E2
Gracey 11/12 Curette
Rigid #9



SG13/14R9E2
Gracey 13/14 Curette
Rigid #9



SG5/6R9E2
Gracey 5/6 Curette
Rigid #9




SG7/8R9E2
Gracey 7/8 Curette
Rigid #9

SUGGESTED INSTRUMENTS (not included in the kit)



RBOE
Boehler Rongeurs
double hinged, 16 cm



H4
#4 Mosquito Hemostat,
Straight 12 cm



DID YOU KNOW?

All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like Instrument Lubricant Spray (**ILS**), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilized in the open position.**

SUGGESTED BLADES AND SUTURES (not included in the Kit)



PSNPTFE1
3-0 PTFE, Undyed Reverse Cutting 19 mm
3/8 circle Reverse Cutting 45 cm



40-815
15 Sterile Scalpel Blades, QTY: 100/box,
Surgical Grade Stainless Steel



MB69
Full Radius Round Scalpel Blades, Sterile,
QTY: 12/box, Surgical Grade Stainless Steel



40-815C
15C Sterile Scalpel Blades, QTY: 100/box,
Surgical Grade Stainless Steel



MB67
Curved Edge Mini Scalpel Blades, Sterile,
QTY: 12/box, Surgical Grade Stainless Steel

Knowledge is power



Subscribe to the **gIDE**
Video Streaming Platform

\$399.00
per Year

\$49.00
per Month



Evidence Based **Video**
Education at its best.



Subscribe here

- 550+ Clinical Videos on Web and App
- New Dental Videos Added Weekly
- Implants, Periodontics, Esthetics & Digital
- For General Dentists & Specialists
- Reduced Pricing on Classroom Training

The gIDE institute is the premiere provider of online dental education with members and online users in over 100 countries and a world renowned faculty of clinicians and researchers from US, Europe and Australasia. gIDE has mastered the technique for delivering high quality digital media - Lecture-On-Demand, Clinical Videos, Expert interviews, Live Broadcasts over the internet to dental professionals, study clubs and societies to specifically meet the educational needs of practicing dentists.



www.gidedental.com

gIDE by the numbers:

108+
Days of
classroom
training

40,000+
Newsletter
subscribers

25,000+
Monthly
visitors

10,000+
Registered
dentists

2,200+
Master clinician
program
graduates

100+
Countries of
members

For more information visit **HuFriedyGroup.com**

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | HuFriedyGroup.com

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC its affiliates or related companies, unless otherwise noted.

©2024 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-943/2407