


# Safe-Flo™ Plus Saliva Ejectors

REF M1002WHWH / M1003CLWH / M1004WHGR / M1005WHPK

## INSTRUCTIONS FOR USE

1. Attach the valve/saliva ejector combination onto the end of the low volume evacuation tube.
2. Bend to the desired position and then place in the mouth.
3. After use, dispose of the valve/ejector combination in the standard trash.



### CROSSTEX



Crosstex International, Inc.  
6789 W. Henrietta Road, Rush, NY 14543 USA  
585.359.0130 | 631.582.6777 | HuFriedyGroup.com

Made in USA

M100P-IFU-EN / REV B - 01/2024

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.

US Pat. No. 8,256,464, 8,714,200, 8,752,579 and foreign counterparts.

©2024 Hu-Friedy Mfg. Co., LLC. All rights reserved.



For professional use only.  
Not made with natural rubber latex.