



DUWL ROUTINE TREATMENT

Liquid Ultra™

IMPORTANT!

**Liquid Ultra™ Solution
Currently in the Waterlines.**

**Flush ALL Dental
Unit Waterlines
per the IFU*
before patient
treatment!**

*Instructions for Use (IFU) can be found on the product box and on the website. For assistance, contact HuFriedyGroup Customer Care 800.483.7433

EPA Reg. 74245-6 and EPA Reg. 74245-7
EPA Est. No. 73516-MD-001

DS503L-HT / REV B - 11/2023

CROSSTEX



DUWL ROUTINE TREATMENT

Liquid Ultra™

IMPORTANT!

**Liquid Ultra™ Solution
Currently in the Waterlines.**

**Flush ALL Dental
Unit Waterlines
per the IFU*
before patient
treatment!**

*Instructions for Use (IFU) can be found on the product box and on the website. For assistance, contact HuFriedyGroup Customer Care 800.483.7433

EPA Reg. 74245-6 and EPA Reg. 74245-7
EPA Est. No. 73516-MD-001

DS503L-HT / REV B - 11/2023

CROSSTEX