

Nebraska Scalers



SN128 Double-End Anterior

The SN128DE includes the Nebraska straight terminal shank with a thinner blade for improved anterior interproximal access. The SN128/Nevi1 Scaler provides the SN128 tip paired with the Nevi1 discoid for excellent anterior stain removal.

- Elongated sickle end for easier access in difficult to reach areas

SN137 Anterior

The SN137 includes a small anterior sickle paired with an anterior curette for improved scaling efficiency.

- A universal H5 sickle on one end for interproximal scaling
- A universal anterior curette on the other end for facial, lingual and cervical areas

Also available in Mini Five™ Design

SN137M Mini Five Design offers a longer shank and shorter, thinner toe on the curette end to provide enhanced access and adaptation in deep or narrow pockets.

SN135 Posterior

The SN135 has a hybrid blade design, offering a rounded back like a curette and a pointed toe of a sickle.

- A rounded back for less trauma when working subgingivally
- A pointed toe for difficult-to-reach areas

Available in EverEdge™ Technology - a proprietary metal finish that stays sharp longer, which means less time spent sharpening.

Nebraska Scalers are available in the following handle options:

Anterior

#128 Double-End | EverEdge™ 2.0

Harmony	SN128DEXE2
#9	SN128DE9E2
#G9, Steel Grip	SN128DEG9E2
#C8, Lavender	SN128DEC8E2

#128 / Nevi1 Discoid | EverEdge™ 2.0

#9	SN128/NEV19E2
----	----------------------

#128 Single End | Immunity Steel™

#8 Gray	SN1288
#7	SN1287
#6	SN1286
#2	SN128

#137 Double-End | EverEdge™ 2.0

Harmony	SN137XE2
#9	SN1379E2
#C8, Lavender	SN137C8E2

Anterior

#137 Double End | Immunity Steel™

#8 Gray	SN1378
#7	SN1377
#6	SN1376

#137 Mini Five™ Double End | Immunity Steel™

#8 Gray	SN137M8
#7	SN137M7
#6	SN137M6

Posterior

#135 | EverEdge™ 2.0

Harmony	SN135XE2
#9	SN1359E2
#G9, Steel Grip	SN135G9E2
#C8, Lavender	SN135C8E2

#135 | Immunity Steel™

#8 Gray	SN1358
#7	SN1357
#6	SN1356

Visit us at **HuFriedyGroup.com**

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.

©2025 Hu-Friedy Mfg. Co., LLC. All rights reserved. HF-573/2510