

# PWR Air

## Air Polisher



Power  
Equipment



# Introducing the **PWR Air**

The PWR Air delivers the power of therapeutic air polishing to every clinician supporting effective supra- and subgingival biofilm disruption, stain removal, and implant cleaning for every patient.

Learn more about PWR Air Polishing at  
[HuFriedyGroup.com/PWR-ON](https://www.hufriedygroup.com/PWR-ON)



PWR Air  
**APAIRPOL**

# Therapeutic **Air** Polishing

## **PWR Air**

Use PWR Air to remove biofilm and tough stains supra- and subgingivally. PWR Air Polishing also helps remove trapped food debris and plaque build-up from between your patients' teeth.

## **What's in the Box**

PWR Air • Air Polishing Device		APAIRPOL
PWR Foot Pedal	<b>UPFTPEDAL</b>	PWR Air K9 Wrench Tip <b>APK9</b>
PWR Air Powder Chamber Prophy	<b>APPCPROPHY</b>	PWR Air Unit Water Supply Hose <b>APH20HOSE</b>
PWR Air Powder Chamber Perio	<b>APPCPERIO</b>	PWR Air Unit Supply Hose <b>APAIRHOSE</b>
PWR Air Powder Chamber Cap	<b>APPCCAP</b>	PWR Air O-Ring Set <b>APORINGS</b>
PWR Air 120° Handpiece	<b>APH120</b>	PWR Air Irrigation Circuit Disinfection <b>APICDKIT</b>
PWR Air Perio Handpiece	<b>APHPPERIO</b>	PWR Air Unit Water Filter <b>APFILTER</b>
PWR Air Perio Subgingival Tips (40/box)	<b>APPERIOTIPS</b>	PWR Air Cleaning Needle Kit <b>APNDLKIT</b>

# PWR Air • Air Polishing

The PWR Air utilizes PWR Prophylaxis Powder combined with air and water to supra- and subgingivally remove biofilm, and the tough stains caused by everyday foods and beverages. Air polishing also helps remove interproximal biofilm and food debris.

## Three Power Modes

Customize your procedures by adjusting the intensity (low, medium or high) according to your patient's sensitivity, type of deposits, and the area being treated—ensuring a targeted, effective, and comfortable approach..

## No Powder Function

Converts the handpiece into an air-water syringe—streamlining rinsing and drying without changing instruments.





Air Polishing is more accessible for biofilm removal in areas of pocketing and around implants. By using the PWR Air, we are able to provide thorough hygiene treatment in areas that are difficult to reach with conventional instrumentation.

**DR. EDWARD P. ALLEN, DDS**



## Slim Stainless Steel Air-Polishing Handpieces

The PWR Air includes three different slim air polishing handpieces: 90°, 120°, and Perio. Being made from steel and resilient to cracking, you can practice with the power peace of mind: fewer interruptions and less time spent replacing handpieces.

### Supra- and Subgingival

The PWR Air allows you to switch easily during treatment between supragingival (Prophy) or subgingival (Perio) biofilm and stain removal with the simple push of a button. The PWR Air automatically performs a flush to clean the internal circuit of any previously used powder, ensuring operational efficiency and ease of use.

### Perio Function

The one-touch Perio Function of the PWR Air reduces the air pressure for a safe and comfortable subgingival treatment.

### Refill Function

The PWR Air Refill Function streamlines maintenance by allowing powder chamber removal and refilling without powering down. This eliminates downtime, optimizes workflow, and enhances operational efficiency for a more seamless practice experience.

### Continuous Flow of Air

PWR Air Device provides a continuous flow of air which helps to prevent clogging. Practice with the power peace of mind: less interruptions and less time spent trying to unclog lines.



### Subgingival Perio Tip

Flexible and soft, the disposable Subgingival Perio Tip (included with the PWR Air) comes sterile and easily attaches to the PWR Perio Handpiece.

### Did You Know?

The PWR Air automatically enters energy-saving mode after 20-minutes of inactivity, reducing power consumption and enhancing device safety.

# PWR Air • Air Polishing Handpieces & Powders

## The PWR of Three

Achieve Best in Practice access with the PWR Air 90° and 120° Angled Handpieces for supra- and subgingival use in periodontal pockets up to 5 mm deep. The PWR Air Perio Handpiece equipped with the disposable sterile PWR Subgingival Perio Tip (included with the PWR Air) provides optical access in pockets deeper than 5 mm. PWR Air Polishing Handpieces easily attach with just a click and are safely secured, thanks to a dedicated security connection system.

### Did You Know?

All PWR Air Polishing Handpieces are resilient to cracking, being made from steel. Practice with the power peace of mind: fewer interruptions and less time spent replacing handpieces.



PWR Air Perio Subgingival Tips  
**APPERIOTIPS**  
QTY: 40



PWR Air 90° Handpiece  
**APHP90**



PWR Air 120° Handpiece  
**APHP120**

Included with the PWR Air



Included with the PWR Air

PWR Air Perio Handpiece  
**APPPERIO**

## PWR Prophylaxis Powders

PWR Prophylaxis Powders remove supra- and subgingival biofilm, stains and pigmentations along with treating peri-implantitis and cleaning of prosthesis and restorative materials. Choose from PWR Perio Powder (Glycine) and PWR Classic Powder (Sodium Bicarbonate), both coming in mint and neutral flavors. Our specially developed PWR Prophylaxis Powders also come with hydrophobic effect, which helps decrease clogging in the PWR Air.



Powder	Flavor	Type	SKU
<b>Perio</b>	Mint	Glycine	<b>AP-PM</b>
<b>Perio</b>	Neutral	Glycine	<b>AP-PN</b>
<b>Classic</b>	Mint	Sodium Bicarbonate	<b>AP-CM</b>
<b>Classic</b>	Neutral	Sodium Bicarbonate	<b>AP-CN</b>

### Tips for PWRful Procedures

PWR Prophylaxis Powders have you covered no matter the situation.

	Classic	Perio
Supragingival Biofilm Removal	●	●
Removal of Stains & Pigmentations on Enamel	●	
Subgingival Biofilm Removal		●
Periodontal Pocket Detoxification		●
Implant Maintenance		●
Peri-Implantitis Treatment		●
Fixed Orthodontic Appliances		●
Cleaning of Prosthesis & Restoration Materials		●
Restorations Maintenance		●

# Ensure Your **PWR Devices** Are Maintained for Optimal Performance



## Our Best in Practice EverCare Service & Repairs **PWR Program.**



For More Membership Information, Visit  
[EverCareService.com](https://EverCareService.com)