

IMS™ Signature Series™ Cassette Clip

Part Code: IM1002
Stainless Steel
Contains: 6 (Six)
Clean and sterilize before use.



www.hufriedygroup.com/Reprocessing/
1-800-Hu-Friedy; info@hu-friedy.com



134°C
Sterilizable in steam sterilizer
at 134°C / 273°F for 4 min



Washer-disinfector for thermal
disinfection possible


All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC
its affiliates or related companies, unless otherwise noted.
©2023 Hu-Friedy Mfg. Co., LLC. All rights reserved.

∴HuFriedyGroup
The Best In Practice

INSTRUCTIONS

Use clip on lid of cassette.

CAUTION: Overextending, bending or attempting to modify
can cause breakage. Clip is recommended for use with IMS™
Signature Series™ Cassettes only.



Hu-Friedy

Hu-Friedy Mfg. Co., LLC
3232 N. Rockwell Street
Chicago, IL 60618 | USA



PLASTIC
BAG

B2939/1023