



REF IMS-1379

QTY 1

Infinity Series™
Barcode Plate
Plaque de code à barres

..HuFriedyGroup
The Best In Practice



Sterilizable in steam sterilizer
at 134 °C / 273 °F for 4 min



Washer-disinfectors for thermal
disinfection possible



www.hufriedygroup.com/Reprocessing/
1-800-Hu-Friedy; info@hu-friedy.com

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC
its affiliates or related companies, unless otherwise noted.
©2025 Hu-Friedy Mfg. Co., LLC. All rights reserved.

Hu-Friedy



Hu-Friedy Mfg. Co., LLC
2501 Sanders Rd.
Northbrook, IL 60062 | USA



N3290 / REV 3 - 2508