

Color Code Instrument Rings
Anneaux de codage couleur pour instruments

QTY 50 | Black / Noir



+D67118275319



Sterilizable in steam sterilizer
at 134 °C / 273 °F for 4 min

All company and product names are trademarks of
Hu-Friedy Mfg. Co., LLC, its affiliates or related
companies, unless otherwise noted.

©2025 Hu-Friedy Mfg. Co., LLC. All rights reserved.

Hu-Friedy



Manufactured for:
Hu-Friedy Mfg. Co., LLC
2501 Sanders Rd.
Northbrook, IL 60062 | USA
HuFriedyGroup.com

IMS-1289

N3062 / REV 3 - 2511

..HuFriedyGroup
The Best In Practice

INK COLORS



BLACK