

**IMS™ Color Code Tabs**  
**Onglets de code couleur**



Sterilizable in steam  
sterilizer at 134 °C /  
273 °F for 4 min

**IMS-12913N**

5 Tabs/Brown | 5 onglets/Brun  
N2888 / REV 3 - 2607



Manufactured for:

**Hu-Friedy Mfg. Co., LLC**

2501 Sanders Rd.

Northbrook, IL 60062 | USA

[HuFriedyGroup.com](http://HuFriedyGroup.com)

**..HuFriedyGroup**



**+D6714681181B**