

Atraumair Forceps Kit

ATRAUMFORCEPS

The Right Pattern Changes the Approach

Six Specialized Forceps. Targeted Extraction Coverage.

Clinicians may rely on a small number of familiar forceps for a wide range of extraction cases. However, comprehensive extraction coverage is not simply about having more instruments. It is about having the appropriate forceps pattern for the clinical application. The **Atraumair Forceps Kit** brings together six complementary standard, apical, root, and cowhorn patterns selected to address common upper and lower extraction needs.

Atraumair Forceps are designed for optimal access and enhanced crown and root engagement. Enhanced beak geometries and longitudinal beak serration provide a firm, secure grasp of the crown and root, reducing the risk of crushing or fracturing the tooth. The unique hole design reduces the weight of the forceps for increased comfort and control, while the matte finish enhances contrast and reduces light reflection.


From upper incisors and roots to lower anterior teeth, premolars and molars, the Atraumair Forceps Kit gives clinicians targeted options across a broad range of extraction cases. Whether standardizing a procedural setup or expanding an existing forceps collection, this thoughtfully assembled kit provides the patterns needed to approach each application with greater control.

Coverage:

FAF35XS	35 Apical Forceps	Upper Canines and Premolars
FX79XS	79 Forceps	Lower Molars
F1XS	1 Standard Forceps	Upper Incisors
FAF36XS	36 Apical Forceps	Lower Premolars, Canines and Incisors
F23XS	23 Cowhorn Forceps	Lower Molars
FX51XS	51 Root Forceps	Upper Roots, Serrated

Visit us online at [HuFriedyGroup.com](https://www.HuFriedyGroup.com)

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.
©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-1323/2609



FAF35XS

FX79XS

F1XS

FAF36XS

F23XS

FX51XS