

3rd Molar Essentials Kit

3RDMOLARESS

Third Molar Readiness Starts with the Right Surgical Essentials


Third molar procedures can require tissue reflection, tooth luxation, root retrieval, extraction, and final bone refinement. The **3rd Molar Essentials Kit** brings together the core instruments needed to support these key stages in one coordinated procedural setup.

Carefully curated for third molar procedures, the kit combines upper and lower apical forceps, miniature right and left Cryer elevators, a slim luxating elevator, a Cogswell elevator, a periosteal elevator, a curved hemostat, and a bone file. Together, these instruments support tissue elevation, controlled tooth and root mobilization, extraction, and bone refinement.

Whether standardizing an existing surgical setup or building a dedicated third molar tray, the 3rd Molar Essentials Kit provides a focused selection of instruments designed to support the key stages of the procedure. By organizing procedure essentials into one solution, the kit helps streamline tray setup and supports consistency across operators.

Kit Contents

BF1X6	1X Miller-Colburn Bone File
E21	Miniature Cryer, Left
E22	Miniature Cryer, Right
FAF150	150 Apical Forceps, Upper Universal
FAF151	151 Apical Forceps, Lower Universal
EB	B Cogswell Elevator
H2	#2 Curved Kelly Hemostat
ELSLIM3CSX	Slim Luxating Elevator, 3mm, Short Curved, Black Line
P96	#9 Molt Periosteal Elevator



Procedural Coverage:

- Tissue Elevation
- Tooth Luxation
- Root Retrieval
- Upper and Lower Extractions
- Tissue Management
- Socket Debridement
- Bone Refinement

Visit us online at [HuFriedyGroup.com](https://www.HuFriedyGroup.com)

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.
©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-1328/2609

 **HuFriedyGroup**
The Best In Practice