



HU-FRIEDY KEY OPINION LEADER

INSTRUMENT KIT LIST



Dr. Bertrand Khayat, D.D.S., MSD

Dr. Khayat is a leading practitioner and global lecturer in Microsurgical Endodontic techniques. He has given over 200 presentations and has published the results of his research in many endodontic journals, including the Journal of Endodontics (JOE). Active in many associations, he maintains his private practice in Paris, France.

His current research interests include new techniques for improved canal preparation. Developed in conjunction with Hu-Friedy, his new hand instrument kit was designed to facilitate these procedures.

DR. BERTRAND KHAYAT

Endodontic Surgical Kit

Part Code: **IMKHAYAT2**

DESCRIPTION	PART CODE
Large Signature Series® 20 Instrument Cassette, Blue*	IM4208
#5 Scalpel Handle, European Style, Round	10-130-05E
Microsurgical Scalpel Blade/Mirror Handle, Stainless Steel	SHDPV
Micro Mirror Round 2.5mm	MM4
NEW – Khayat Apical Plugger with Satin Steel® Handle	PLGBK6
NEW – Khayat Filling Instrument with Satin Steel Handle	PF18ABK6
85 Lucas Surgical Curette, Spoon Shape, Short Satin Steel Handle	CL856
Velvart Apical Plugger with Satin Steel Handle	PLGVEL56
NEW – Perma Sharp® General Tissue Pliers, 15cm/6"	TP5072
PR-3 Prichard Periosteal, with Satin Steel Handle	PPR3S6
13/14 Columbia Curette, Universal with Satin Steel Handle	SC13/146
34/35 Jacquette Scaler, Miniature, Posterior, with Satin Steel Handle	SJ34/356
89/92 Carver Discoid-Cleoids, with Satin Steel Handle	CD89/926
DG 16/17 Endodontic Explorer, with Satin Steel Handle	EXDG16/176
Cawood-Minnesota Retractor	CRM2
NEW – Perma Sharp Needle Holder, Castroviejo Curved, 14cm/5", Round Handle	NH5021RC
14 LaGrange Double-Curved Super Cut Scissors	S14SC

All instruments sold in full kit sets are also available for individual purchase.

**Kit includes hinged instrument clips.*



How the best perform

Hu-Friedy

NEW INSTRUMENTS



KHAYAT APICAL PLUGGER

The Khayat Apical Plugger, with its 4.5mm and 6.0mm working ends, is ideal for packing retrofilling materials. It was specifically designed for use in deep canal preparation.


PART CODE: PLGBK6



KHAYAT FILLING INSTRUMENT

The Khayat Filling Instrument has one angled, and one straight spatula end. The straight spatula end is excellent for hard to reach areas in complex surgical procedures.

PART CODE: PF18ABK6



PERMA SHARP GENERAL TISSUE PLIERS

The Perma Sharp Tissue Pliers combine the advantages of "plain" carbide inserts in limiting soft tissue trauma with the ergonomic benefits of a sculpted handle.

PART CODE: TP5072



PERMA SHARP CASTROVIEJO NEEDLE HOLDER

The Perma Sharp Castroviejo Needle Holder combines a curved beak with the maneuverability provided by its round handle. (Recommended for suture sizes 5.0, 6.0, 7.0.)

PART CODE: NH5021RC

Visit us online at [Hu-Friedy.com](https://www.hu-friedy.com)

©2026 Hu-Friedy Mfg. Co., LLC. All rights reserved. M3044/0617 [630]

Hu-Friedy