



HU-FRIEDY KEY OPINION LEADER

PROF. GIULIO RASPERINI
PERIODONTAL SURGERY SET + TUNNELLING SET

Prof. Giulio Rasperini received his degree in dentistry and specialised in Orthodontics.



He is an active member of the Italian Society of Periodontology, the European Academy of Esthetic Dentistry, and is an ITI fellow.

He serves on the Editorial Boards for several publications, including the International Journal of Periodontics and Restorative Dentistry and the Journal of Implant and Advanced Clinical Dentistry.

Dr. Rasperini has written several publications focused on Periodontology and Implantology, and has won awards for his research—most recently, the Earl Robinson Periodontal Regeneration Award from the American Academy of Periodontology.

He is a professor in Milan, Italy and Michigan, United States, and maintains a private practice specialising in periodontics and implant therapy.

“ The evolution of the surgical techniques in Periodontology in the last years leads as consequence, to an evolution of the surgical instruments. The introduction of Ultrasonic or Piezo tips for osseo-resective procedure, loops and operative microscope redesign the surgical armamentarium. This Kit has been thought for all the most recent periodontal procedures. ”


– PROF. GIULIO RASPERINI

Cassette Color Coding

- IMEDIN180S = Red
- IMEDIN181S = Grey
- IMEDIN182S = White
- IMEDIN183S = Orange
- IMEDIN184S = Purple (new)
- IMEDIN185S = Yellow
- IMEDIN186S = Lavender
- IMEDIN187S = Black (new)
- IMEDIN188S = Blue
- IMEDIN189S = Green
- IMEDIN1811S = Ocean Blue (new)

IMEDIN188S

Next Generation IMS Cassette Signature DIN 18 Instruments, blue. With accessory area and clips for hinged instruments.



Hu-Friedy IMS Instrument Management System

By standardising your instruments with IMS, additional time and revenue is generated while extending the useful life of your Hu-Friedy instruments. Dental professionals state that by implementing IMS, they save 5 to 10 minutes per procedure allowing the staff to spend more time on value-added activities. In addition, safer transportation of soiled instruments reduces risk of injury and exposure to blood borne pathogens.



MIR5DSHDE
 High Definition CS Mirror #5
 double-sided, European Stem,
 1pc/pack

WHY OUR HD MIRRORS ARE SUPERIOR

- **113% reflection factor** for exceptional image clarity
- **38,5% brighter** than rhodium coated mirrors
- **50% brighter** than other front surface mirrors
- **Scratch-resistant** for product longevity
- **Large ergonomic handle** options



MHE6
 Mirror Handle hdl #6
 European Thread



DP17
 Tweezers #17


XP23/UNC6
 Expro North Carolina #23 hdl #6
 15-14-13....3-2-1



PCORTLINI
 Probe Cortellini hdl #30
 15-....0,5

COLORVUE BIOTYPE PROBES

- Colorvue Biotype Probe tips are designed for up to 30 uses
- Colorvue Biotype Probes fit in Hu-Friedy's ergonomic, Satin Steel handle that can be reused or recycled
- The resin tips are safe and gentle to use on all patients
- Each colored tip is also numbered to aid in differentiation between colors



PBTPKIT12
 Biotype Probe Starter kit
 (3 handles + 12 tips)



PQ2N6
 Bifurcation Probe Nabers #2N hdl #6
 Quilix 3-6-9-12



10-130-05E
 Scalpel Handle #5E Europ.



SHDPV
 Scalpel Handle #5E Europ.



TPSLCOSM
 Tissue Pliers Micro Cooley hdl
 #SinusLine straight modified 18cm

PROF. GIULIO RASPERINI
KIT CHIRURGIA PARODONTALE



PP5590SIM
 Periosteal Simion hdl #6



P24GSIC6
 Periosteal #24 Sicilia #6 long hdl #6



PH26M
 Periosteal micro surgery #2MBH hdl #6

EVEREDGE® TECHNOLOGY

Sharpen Less Often

Revolutionary new stainless steel alloy is super-durable. Stays sharper longer. That means less time sharpening and less hand fatigue.

Proprietary Process

Proprietary heat treat and cryogenic processing ensure that the superior edge retention and wear characteristics of EVEREDGE® Technology will last the entire life of the instrument. It's not a superficial coating - EverEdge Technology® scalers can be sharpened again and again for your best instrument value.

Diamond Knurl

Unique diamond knurl pattern with large diameter handle for a confident grasp. Optimal weight for reduced hand fatigue.



SAS7/897
 Curette Gracey #7/8 hdl #9 MF oral/buccal, green



SAS11/1293
 Curette Gracey #11/12 hdl #9 MF mesial, orange



SAS13/1498
 Curette Gracey #13/14 hdl #9 blue MF distal




SGF46
 Curette Goldman-Fox #4 hdl #6 universal



EXC31W
 Excavator Endo #31W hdl #41




C36/376
 Back Action Bone Chisel #36/37 hdl #6



NHM-5026R
 Needle Holder Castroviejo, round diamond dusted beak, straight, 18cm



SCORT
 Scissors Cortellini straight 12cm



IME14DIN5_
Next Generation IMS Cassette
Signature 1/4 DIN, 5 Instruments

Cassette Color Coding

- IME14DIN50 = Red
- IME14DIN51 = Grey
- IME14DIN52 = White
- IME14DIN53 = Orange
- IME14DIN54 = Purple (new)
- IME14DIN55 = Yellow
- IME14DIN56 = Lavender
- IME14DIN57 = Black (new)
- IME14DIN58 = Blue
- IME14DIN59 = Green
- IME14DIN511 = Ocean Blue (new)

The Black Line Collection of over 50 instruments features:

- Performance engineered coating for a harder, smoother surface for **optimal edge retention and enhanced lubricity**
- Distinct black finish for **enriched contrast and visual acuity at the surgical site and underlying tissue**
- **Reduced light reflection** afforded by a matte finished handle and black working ends
- Unique, smooth, large diameter lightweight handle for **maximum comfort, reduced hand fatigue and increased control**



TKN1SX
Tunneling Knife #1 Black Line, Short Handle



TKN2SX
Tunneling Knife #2 Black Line, Short Handle



PPAELSX
Periosteal Elevator Allen, Black Line, Short Handle



FQUINCEY1X
Tunneling Instrument GDQ #1, Black Line



FQUINCEY2X
Tunneling Instrument GDQ #2, Black Line

