



HU-FRIEDY KEY OPINION LEADER

DR. GIORGIO TABANELLA
PERIOSURGICAL KIT

Dr. Giorgio Tabanella

Dr. Tabanella is a Diplomate of the American Board of Periodontology. He graduated from the University of Southern California-Los Angeles-USA where he obtained the Certificate in Periodontics as well as the Master of Science in Craniofacial Biology.

He was awarded "Outstanding Periodontal Graduate Student Research In Surgery And Implantology" by the California Society of Periodontists for his study on bone remodeling.

He is the Scientific Coordinator of the International Implant Conference, and reviewer and author of original articles published in peer-reviewed journals such as the International Journal of Oral and Maxillofacial Implants, Journal of American Dental Association, Journal of Inflammation, World Journal of Stomatology, Clinical Implant Dentistry and Related Research.

Dr. Tabanella lectures in Europe as well as in the US on implant and periodontal surgery, aesthetics, peri-implant bone remodeling, soft and hard tissue reconstruction.

His research focuses on tissue regeneration and augmentation as well as peri-implant bone loss. He maintains a private practice in Rome, Italy.

www.tabanellaimplantologia.it*



TABANELLA1
Tabanella Micro Elevator

Sharp micro elevator for rapid and minimally invasive flap and papilla elevation. Can also be used for removal of secondary flaps.



TABANELLA2
Tabanella Universal Bone File

Sharp Universal bone file designed to dissect thick fibrous tissue, harvest small bony chips (ideal for regenerative surgery) and remove ankylosed root tips.



TABANELLA3
Tabanella Bone Chisel 1

Sharp periodontal bone chisel designed for interproximal osteoplasty during osseous surgery.



TABANELLA4
Tabanella Bone Chisel 2

Sharp periodontal bone chisel designed to perform interproximal ostectomy during osseous surgery on the mesial and distal aspect of the root surface.



MIR5DSHD
HD Mirror #5, Double Sided,
US-Thread



MH6
Mirror Handle hdl #6 US-Thread



DP18
Tweezers #18



CRM
Cheek Retractor Minnesota



XP23/UNC6
Expro- N Carolina #23 hdl #6
15-14-13....3-2-1

EVEREDGE® TECHNOLOGY

A revolutionary stainless steel alloy provides optimal cutting edge.

Proprietary Process

Proprietary heat treat and cryogenic processing ensure that the superior edge retention and wear characteristics of EVEREDGE® Technology will last the entire life of the instrument.

Diamond Knurl

Unique diamond knurl pattern with large diameter handle for a confident grasp. Optimal weight for reduced hand fatigue.



S13K/TG9
Curette/Chisel#13K/TG DE hdl #9



SPR1/29
Curette #PR1/2 Prichard hdl #9




SC4R/4L9
Curette Columbia #4R/4L hdl #9



10-130-05E
Scalpel Handle #5E European Style



10-256-15
Scalpel Blades #15C, 100 pcs/pkt




TPSLCOSM
Tissue Pliers Micro Cooley hdl #Sinus-
Line straight 18cm



P20
Periosteal Hirschfield #20



S5081
Scissors Goldman-Fox #5081 angled,
PermaSharp 12,5cm



The Black Line Collection of over 50 instruments features:

- Performance engineered coating for a harder, smoother surface for **optimal edge retention and enhanced lubricity**
- Distinct black finish for **enriched contrast and visual acuity at the surgical site and underlying tissue**
- **Reduced light reflection** afforded by a matte finished handle and black working ends
- Unique, smooth, large diameter lightweight handle for **maximum comfort, reduced hand fatigue and increased control**



PPBUSERX
Periosteal Buser Black Line



PPBUSER12X
PJO/1+PJO/2 combined, 13KL, Black Line



PPR3X
Periosteal Prichard PR-3, Black Line



HKRC
Hemostat Kelly,
curved 16cm



FX36E
Forceps #36
European IM
premolar



NH5024R
Needle Holder Castroviejo
#5024R round PermaSharp





RBEY
Rongeur Bayer 18cm

PSN493V
Surgical Suture PGA UNDY #5-0A/C-3 P.
reverse cut 3/8crl

PSN8683P
Surgical Suture Polypropylene #4-0NA/C-6
reverse cut 3/8crl

PSN8698P
Surgical Suture
Polypropylene
#5-0NA/C-3 P.
reverse cut 3/8crl



TRE120M
Trephine Bur
inner diameter,
right angle 12mm

HU-FRIEDY NEW GENERATION IMS CASSETTES

Hu-Friedy's Instrument Management System (IMS) is designed to keep instruments organised and intact from cleaning to chairside, providing opportunity for:



- 5-10 minutes saved per procedure
- Standardised procedures
- Increased staff safety
- Increased instrument life
- Focus on revenue-generating activities
- Reduction in overtime paid
- Potential to see 1-2 additional patients per day

IMEDIN188S

Next Generation IMS Cassette
Signature DIN 18 Instruments, blue.
With accessory area and clips for
hinged instruments.

Cassette Color Coding

- IMEDIN180S = Red
- IMEDIN181S = Grey
- IMEDIN182S = White
- IMEDIN183S = Orange
- IMEDIN184S = Purple (new)
- IMEDIN185S = Yellow
- IMEDIN186S = Lavender
- IMEDIN187S = Black (new)
- IMEDIN188S = Blue
- IMEDIN189S = Green
- IMEDIN1811S = Ocean Blue (new)




PERIODONTAL



RESTORATIVE



ORTHODONTIC




INSTRUMENT
MANAGEMENT



ENDODONTIC



SURGICAL



ULTRASONIC
SCALING



Please visit us at www.Hu-Friedy.eu

How the best perform

