

BEAUTIFUL BIOMETRIC RESTORATIONS

| DR. CLARENCE TAM



FIGURE 1

PREOPERATIVE VIEW

Defective amalgam and composite patch restorations which need to be replaced.



FIGURE 2

Amalgams removed and marginal ridge composite base layer shaped in the distal prep using TNTAM2.



FIGURE 3

Brown tint delicately applied with the ultra-pointed fine end of TNTAM1.



FIGURE 4

Final layers of composite placed and sculpted with the conical end of TNTAM1.



FIGURE 5

Access for matrix burnishing and initial layers of packable composite made simple using the extended shank and burnisher end of TNTAM3.



FIGURE 6

Composite formed and finessed with the thin paddle of TNTAM2. Final layer of composite placed and shaped with TNTAM1.



FIGURE 7

POSTOPERATIVE VIEW

Beautifully and biomimetically-completed restorations.



FIGURE 8

How the best perform

Hu-Friedy



Dr. Clarence Tam

Dr. Clarence Tam is originally from Toronto, Canada, where she completed her Doctor of Dental Surgery and General Practice Residency at the University of Western Ontario and the University of Toronto, respectively. Clarence's practice is limited to cosmetic and restorative dentistry. She is well-published to both the local and international dental press, writing articles, reviewing and developing prototype products and techniques in clinical dentistry. She frequently and continually lectures internationally.

Clarence is the Chairperson and Director of the New Zealand Academy of Cosmetic Dentistry. She is currently one of two individuals in Australasia to hold Board-Certified Accredited Member Status with the American Academy of Cosmetic Dentistry. Clarence is an Opinion Leader for multinational dental companies Kuraray Noritake, J Morita Corp, Henry Schein NZ, Ivoclar Vivadent, DentsplySirona, 3M ESPE, Kerr, GC Australasia, SDI, Coltene, Triodent and is the only Voco Fellow in Australia and New Zealand. She holds Fellowship status with the International Academy for DentoFacial Esthetics. She is a passionate and approachable individual, and is committed to having an interactive approach with patients in all of her cases to maximize predictability.

DESCRIPTION	PART CODE
Cone Socket Mirror Handle, Satin Steel	MH6
#5 HD Mirror	MIR5HD
Posterior Occlusal Shaper	TNTAM1
Universal Thin PFI	TNTAM2
Curved Ball Burnisher	TNTAM3
Posterior Sickle Scaler, EverEdge 2.0	S204S9E2
Complete Kit	TNTAMKIT*

*Kit contains instruments in 5-instrument cassette



Universal Thin PFI
| **TNTAM2**



Posterior Occlusal Shaper
| **TNTAM1**



Curved Ball Burnisher
| **TNTAM3**



Learn more about StyleSkulpt and Dr. Tam at:
[Hu-Friedy.com/Community/Clarence-Tam](https://www.hu-friedy.com/Community/Clarence-Tam)

©2018 Hu-Friedy Mfg. Co., LLC. All rights reserved. [910] HFL-046/0718

