

Wrap Size Recommendations

| | 10" x 10" | 15" x 15" | 20" x 20" | 24" x 24" |
|---|-----------|-----------|-----------|-----------|
| Signature Series® Cassettes | | | | |
| Exam, Double-Decker™ 8-Instrument, Ultrasonic Insert, Symmetry IQ® & Handpiece Cassette | X | X | | |
| Space Saver Seven™ Cassette | X | X | | |
| Double-Decker 14-Instrument Cassette Δ | | X | X | |
| Small 8 or 10 Instrument Cassette Δ | | X | X | |
| Signa-Stat or Double-Decker 20-Instrument Cassette | | | X | |
| Large 16 or 20 Instrument Cassette | | | X | |
| Oral Surgery Cassette | | | | X |
| Orthodontic Cassette Box | | | | X |
| Large Thin Orthodontic Cassette | | | X | |
| Small Thin Orthodontic Cassette | | | X | |
| Double-Decker Large Ortho Cassette | | | X | |

Δ The Double-Decker 14-instrument cassette and Small cassettes are best with the 15" x 15" sterilization wrap. However, if an office is trying to purchase only one size, the 20" x 20" size is recommended.