




ADA CDT CODING COLLECTION SERIES

Hu-Friedy has created a collection of procedural instrument kits that follow **American Dental Association (ADA) – Current Dental Terminology (CDT) Coding**. Start pairing patients and insurance codes more effectively to move from standardized care to personalized care.

Begin each preventative appointment with an Evaluation Kit and add the supplemental treatment kit based on the patient assessment.


Each kit, uniquely color-coded, offers the highest quality, technology, and ergonomic instrumentation within a small, space-efficient cassette for ensured equipment compatibility and safety. Plastic Test Sticks are included to test for sharpness prior to utilizing instrumentation to provide ergonomic, efficient and quality care to your patients. ADA CDT Coding kits can be ordered easily using a single part code at a significantly reduced bundled cost.

NOTE: Collection Series cassettes can be laser marked with a part code and/or kit description for easy identification. Contact your Hu-Friedy Account Manager for more information.



DID YOU KNOW?

The Mesial 15/16 Gracey evolved from the 11/12, created while clinicians practiced stand up dentistry and the occlusal plane was parallel to the floor. The 15/16 has an extra bend in the functional shank to maintain clinicians' wrists in a neutral position.



Evaluation Kit

KIT CONTENTS

Exam 5 Instrument Cassette, Blue	IM6058
Mirror Handle, SE Satin Steel #6 Handle	MH6
5 Mirror Head, Double Sided High Definition	MIR5DSHD
#5 Explorer, DE Satin Steel #6 Handle	EXD56
UNC Probe, SE Satin Steel #6 Handle	PCPUNC126
#11/12 ODU DE After 5 Explorer, DE Satin Steel #6 Handle	EXD11/12A6
Nabors Probe 3-6-9-12, DE Satin Steel #6 Handle	PQ2N6


| **IM-D0120/80**

Prophylaxis Kit

KIT CONTENTS

Exam 5 Instrument Cassette, Green	IM6059
Plastic Test Stick – Individual	PTS-IND
Nevi 2 Sickle Scaler, DE EverEdge 2.0 #9 Handle	SCNEVI29E2
Nevi 4 Sickle Scaler, DE EverEdge 2.0 #9 Handle	SCNEVI49E2
5/6 Barnhart Universal Curette, DE EverEdge 2.0 #9 Handle	SBH5/69E2
15/16 Gracey Curette, DE EverEdge 2.0 #9 Handle	SG15/169E2
13/14 Gracey Curette, DE EverEdge 2.0 #9 Handle	SG13/149E2

| **IM-D1110**



Pedo Prophy Kit

KIT CONTENTS

Exam 5 Instrument Cassette, Lavender	IM6056
Plastic Test Stick – Individual	PTS-IND
Mirror Handle, SE Satin Steel #6 Handle	MH6
5 Mirror Head, Double Sided High Definition	MIR5DSHD
23/5 Explorer, DE Satin Steel #6 Handle	EXD56
13/14 Columbia Universal Curette, DE EverEdge 2.0 #9 Handle	SC13/149E2
Nevi 3 Sickle Scaler, DE Ever Edge 2.0 #9 Handle	SCNEVI39E2
Nevi 1 Sickle Scaler, DE Ever Edge 2.0 #9 Handle	SCNEVI19E2

| **IM-D1120**



How the best perform



Scaling and Root Planing Kit

KIT CONTENTS

Double Decker, 8 Instrument Cassette, Red	IM9080
Plastic Test Stick – Individual	PTS-IND
3/4 Langer After 5 Universal Curette, DE EverEdge 2.0 #9 Handle	SL3/4AF9E2
Nevi 2 Sickle Scaler, DE EverEdge 2.0 #9 Handle	SCNEVI29E2
Nevi 4 Sickle Scaler, DE EverEdge 2.0 #9 Handle	SCNEVI49E2
1/2 Mini 5 Gracey Curette, DE EverEdge 2.0 #9 Handle	SASI/29E2
7/8 After 5 Gracey Curette, DE EverEdge 2.0 #9 Handle	SRPG7/89E2
15/16 After 5 Gracey Curette, DE EverEdge 2.0 #9 Handle	SRPG15/169E2
13/14 After 5 Gracey Curette, DE EverEdge 2.0 #9 Handle	SRPG13/149E2
17/18 Gracey Curette, DE Ever Edge 2.0 #9 Handle	SG17/189E2



| IM-D4341/42

Presence of Gingival Inflammation Kit

KIT CONTENTS

Exam 5 Instrument Cassette, Yellow	IM6055
Plastic Test Stick – Individual	PTS-IND
Nevi 3 Sickle Scaler, DE EverEdge 2.0 #9 Handle	SCNEVI39E2
Nevi 4 Sickle Scaler, DE EverEdge 2.0 #9 Handle	SCNEVI49E2
1/2 Mini 5 Gracey Curette, DE EverEdge 2.0 #9 Handle	SASI/29E2
15/16 After 5 Gracey Curette, DE EverEdge 2.0 #9 Handle	SRPG15/169E2
13/14 After 5 Gracey Curette, DE EverEdge 2.0 #9 Handle	SRPG13/149E2



| IM-D4346

Debridement for Comprehensive Evaluation Kit

KIT CONTENTS

Exam 5 Instruments Cassette, Orange	IM6053
Plastic Test Stick – Individual	PTS-IND
Nevi 4 Sickle Scaler, DE EverEdge 2.0 #9 Handle	SCNEVI49E2
17/18 McCalls Universal Curette, DE EverEdge 2.0 #9 HdI	SM17/189E2
1/2 Rigid Mini 5 Gracey Curette, DE EverEdge 2.0 #9 Handle	SASI/2R9E2
15/16 Rigid After 5 Gracey Curette, DE EverEdge 2.0 #9 Handle	SRP15/16R9E2
13/14 Rigid After 5 Gracey Curette, DE EverEdge 2.0 #9 Handle	SRP13/14R9E2




| IM-D4355

DID YOU KNOW?

Anna Pattison Curettes are modified from the original Gracey Curettes. They have shorter, thinner, smaller blades honed at a 60 degree angle with a modified rigid shank specially designed to allow easier access and easier insertion under tight tissue.

Anna Pattison Curette

Gracey Curette



Periodontal Maintenance Kit

KIT CONTENTS

Double Decker, 8 Instrument Cassette, Gray	IM9081
Plastic Test Stick – Individual	PTS-IND
3/4 Langer Universal Curette, DE EverEdge 2.0 #9 Handle	SL3/49E2
Nevi 2 Sickle Scaler, DE EverEdge 2.0 #9 Handle	SCNEVI29E2
Nevi 4 Sickle Scaler, DE EverEdge 2.0 #9 Handle	SCNEVI49E2
1/2 Pattison Gracey LITE Curette, DE EverEdge 2.0 #9 Handle	SG1/29APE2
7/8 Pattison Gracey LITE Curette, DE EverEdge 2.0 #9 Handle	SG7/89APE2
11/12 Pattison Gracey LITE Curette, DE EverEdge 2.0 #9 Handle	SG11/129APE2
13/14 Pattison Gracey LITE Curette, DE EverEdge 2.0 #9 Handle	SG13/149APE2
#15/16 Gracey Curette, EverEdge 2.0 #9 Handle	SG15/169E2

| IM-D4910

Implant Maintenance Kit

KIT CONTENTS

Exam 5 Instrument Cassette, White	IM6052
Plastic Test Stick – Individual	PTS-IND
Colorvue Probe, UNC 12, SE Satin Steel Handle	PCVPROMONC
Titanium Implant Scaler, 204SD Sickle Scaler	IMP204SDT
Titanium Implant Scaler, Langer 1/2 Universal Curette	IMPLG1/2T
Titanium Implant Scaler, 1/2 Gracey Mini Five Curette	IMPM1/2T
Implacare II Implant Scaler, H6/7 Sickle Scaler	IMPLH6/76

| IM-D6080

