



Akro-Flex™


Create Beautiful Restorations

Akro-Flex™ is an extremely flexible composite instrument which allows clinicians to effortlessly manipulate composite materials during aesthetic restorations. The inclusion of a ductile material allows Akro-Flex to perform as a solid brush.

HYPER-THIN PROFILE

Ultra thin working ends reach narrow interproximal spaces with ease.

The flexible, versatile working ends allow for better visibility as compared to traditional composite instruments.



UNPARALLELED FLEXIBILITY

By incorporating Nickel Titanium, a material known for its super elasticity, Akro-Flex acts as a solid brush.

The resilient working ends are excellent when creating fine anatomical detail with delicate, artistic strokes.

ERGONOMIC HANDLE

The smooth, lightweight handle offers increased control due to the large diameter. It creates an ergonomically friendly option that provides maximum comfort and helps reduce hand fatigue. Reducing hand fatigue can increase the longevity of a clinician's career.

To learn more, visit Hu-Friedy.com/CompositeInstruments

How the best perform

