

**Product Name:** ENZYMAX Concentrated Ultrasonic Detergent and Presoak

**Product Code:** IMS-1222

| <b>Ingredients</b>    | <b>CAS No.</b> | <b>Function</b> |
|-----------------------|----------------|-----------------|
| Water                 | 7732-18-5      | Diluent         |
| Alkanoate salt        | 764-71-6       | Cleaning agent  |
| Glycerin              | 56-81-5        | Stabilizer      |
| Protease enzyme       | 9014-01-1      | Cleaning agent  |
| Amylase enzyme        | 9000-90-2      | Cleaning agent  |
| Acid violet 43        | 4430-18-6      | Colorant        |
| Methylisothiazolinone | 2682-20-4      | Preservative    |
| Benzisothiazolinone   | 2634-33-5      | Preservative    |
| Lime oil              | 8008-26-2      | Fragrance       |
| Lime oil terpenes     | 68917-71-5     | Fragrance       |
| Geranyl acetate       | 105-87-3       | Fragrance       |
| Ethanol               | 64-17-5        | Fragrance       |
| Linalool              | 78-70-6        | Fragrance       |
| Orange terpenes       | 68647-72-3     | Fragrance       |
| Lemon oil             | 8008-56-8      | Fragrance       |
| Geraniol              | 106-24-1       | Fragrance       |
| Octoxynol             | 9036-19-5      | Fragrance       |