

Product Name: Enzymax Earth

Product Code: IMS1332, IMS1332SP, IMS1333

Ingredients	CAS No.	Function
Sodium sesquicarbonate	533-96-0	Diluent
Tetrasodium iminodisuccinate	144538-83-0	Chelating agent
Sodium C10-16 alkylbenzenesulfonate	25155-30-0	Cleaning agent
Sodium metasilicate	6834-92-0	Corrosion inhibitor
Protease enzyme	9014-01-1	Cleaning agent
Cellulase enzyme	9012-54-8	Cleaning agent
d-Limonene	5989-27-5	Fragrance