



# HU-FRIEDY SURGICAL KIT

## CRESTAL SINUS LIFT KIT



PART CODE	DESCRIPTION
<b>IM4162</b>	16 Instrument Signature Series® Cassette, White
<b>10-130-03</b>	3 Scalpel Handle
<b>P96</b>	9 Molt Periosteal Elevator
<b>CM2/4S6</b>	2/4 Molt DE Curette, #6 Handle
<b>TP43</b>	43 Adson-Brown Tissue Pliers
<b>CSEUWES6</b>	Euwe Bone Scoop/Plugger, #6 Handle
<b>PCP65006</b>	20mm Probe 8-10-12-14-16-18-20, #6 Handle
<b>OSTPAL1</b>	1 Palti Osteotome 2.7mm
<b>OSTPAL2</b>	2 Palti Osteotome 3.2mm
<b>OSTPAL3</b>	3 Palti Osteotome 3.7mm
<b>OSTPAL4</b>	4 Palti Osteotome 4.2mm
<b>PFIWDS1</b>	1 Woodson DE Plastic Filling Instrument
<b>PPR3S6</b>	Prichard Periosteal, #6 Handle
<b>IMCRESTKIT</b>	Complete Crestal Sinus Lift Kit

# CRESTAL SINUS LIFT KIT



# CRESTAL SINUS LIFT KIT



**A** **IMN4162**  
**16 Instrument Signature Series® Cassette, White**

IMS® cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.



**B** **10-130-03**  
**3 Scalpel Handle**

Used for cutting gingival tissue and making surgical incisions. Handle is 4.9" long and 13mm wide at base. Compatible with standard scalpel blade sizes.




**C** **P96**  
**9 Molt Periosteal Elevator**

Designed for reflecting and retracting the mucoperiosteum after incisions of the gingival tissue. 4.3mm/8.4mm



**D** **CM2/4S6**  
**2/4 Molt DE Curette, #6 Handle**

Designed for curettage, cyst removal and tooth socket debridement with 3.6mm and 7mm working ends.



**E** **TP43**  
**43 Adson-Brown Tissue Pliers**

Used to grasp and stabilize soft tissue during suturing. 12mm



**F** **CSEUWES6**  
**Euwe Bone Scoop/Plugger, #6 Handle**

Designed for easy transport and non-stick compaction of bone graft material. 3.0mm/7.0mm.



**G** **PCP65006**  
**20mm Probe 8-10-12-14-16-18-20, #6 Handle**

A double-ended instrument with an explorer on one end and a probe on the other, excellent for dental examinations and charting.



**H** **OSTPAL1**  
**1 Palti Osteotome 2.7mm**

Used for preparing bone. 7-10-13-15-18mm



**I** **OSTPAL2**  
**2 Palti Osteotome 3.2mm**

Used for preparing bone. 7-10-13-15-18mm



**J** **OSTPAL3**  
**3 Palti Osteotome 3.7mm**

Used for preparing bone. 7-10-13-15-18mm



**K** **OSTPAL4**  
**4 Palti Osteotome 4.2mm**

Used for preparing bone. 7-10-13-15-18mm



**L** **PFIWDS1**  
**1 Woodson Plastic Filling Instrument**

Designed to reflect and retract the mucoperiosteum after incisions of the gingival tissue with 2.2mm and 3.8mm width.



**M** **PPR3S6**  
**Prichard Periosteal, #6 Handle**

Secure and maintain open position of scissors, hemostats, and other hinged instruments in lid of cassette.