



# HU-FRIEDY SURGICAL KIT

## MICROSURGICAL KIT



PART CODE	DESCRIPTION
<b>10-130-70</b>	7 Mini Scalpel Blade Holder
<b>S13K/13KL6</b>	13K/13KL DE Surgical Chisel
<b>C36/376</b>	36/37 Back-Action DE Chisel
<b>8-905DD</b>	Precision Dressing Forceps, Diamond Dusted
<b>NH5024CORT</b>	Cortellini Micro Castro Needle Holder
<b>PPBUSER6</b>	Buser Periosteal, #6 Handle
<b>SPV</b>	Microsurgical Scissors 18cm
<b>IM8120</b>	12 Instrument Signature Series® Cassette, Red
<b>PCPUNC156</b>	CC SE UNC Probe, 1-15, #6 Handle
<b>MIR5</b>	5 FS CS Mouth Mirror
<b>MH6</b>	Cone Socket Mirror Handle, #6 Handle
<b>IMMICSUR</b>	Complete Microsurgical Kit

# MICROSURGICAL KIT



## MICROSURGICAL KIT



- A 10-130-70**  
**7 Mini Scalpel Blade Holder**  
Designed to hold micro scalpel blades for minimally invasive procedures and microsurgery



- B S13K/13KL6**  
**13K/13KL DE Surgical Chisel**  
Used to recontour periodontal bone or to remove remaining tooth fibers during periodontal surgery with curved working ends 3.2mm and 4.9mm wide.



- C C36/376**  
**36/37 Back-Action DE Chisel**  
Designed for removing bone adjacent to the tooth without causing trauma and is especially useful on the distal of last molars with working ends 3.8mm wide.



- D 8-905DD**  
**Precision Dressing Forceps, Diamond Dusted**  
Used to grasp and stabilize soft tissue, 7" length.



- E NH5024CORT**  
**Cortellini Micro Castro Needle Holder**  
Used to grasp and guide the needle when suturing.



- F PPBUSER6**  
**Buser Periosteal, Satin #6 Handle**  
Designed for reflecting and retracting the mucoperiosteum after incisions of the gingival tissue. 3.8mm/4.3mm



- G SPV**  
**Microsurgical Scissors**  
Microsurgical scissors unique design features allow extremely smooth cutting of fine and coarse tissue with reduced trauma. 18cm



- H IM8120**  
**12 Instrument Signature Series® Cassette, Red**  
IMS® cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.



- I PCPUNC156**  
**CC SE UNC Probe, 1-15, #6 Handle**  
Single-ended, color-coded probe. ergonomic, Satin Steel™ handle. University of North Carolina.



- J MIR5**  
**5 FS CS Mouth Mirror**  
Single-sided Cone Socket mirror. 15/16 in. (24mm) diameter.



- K MH6**  
**Cone Socket Mirror Handle, #6 Handle**  
Used with Single and Double-Sided Cone Socket Mirrors. The blunt end of the mirror handle is ideal for testing sensitivity to percussion and assessing tooth mobility.