


HU-FRIEDY SURGICAL KIT

KEYS SOCKET GRAFTING KIT




PART CODE	DESCRIPTION
IMN4164	16 Instrument Infinity Series™ Cassette, Purple
XP23/86	23/CP-8 Expro
H2	2 Kelly Hemostat, Curved
PLGGR1X	Grisdale Bone Condenser 8-10-13-15-18, Black Line
CM2/4SX	2/4 Molt DE Curette, Black Line
10-130-05	5 Scalpel Handle, Round, Straight
PPR3X	PR-3 Prichard Periosteal, Black Line
NH5037	Halsey Perma Sharp® Needle Holder, Serrated
S5055	Metzenbaum Perma Sharp® Scissors, Curved/Blunt
TP5090	DeBakey Perma Sharp® Tissue Pliers, 6"
DPI	1 Dressing Pliers, Perry
PMPIX	Membrane Placement Instrument, Black Line
MCUPE	Measuring Cup 25 mL w/ Marking
SYBG	Syringe Bone Grafting 7mm
CL856	85 Lucas DE Surgical Curette
PJO/16	S13K/TG Jovanovic 4mm Chisel
IMNKEYSGKT	Complete Keys Socket Grafting Kit

KEYS SOCKET GRAFTING KIT



KEYS SOCKET GRAFTING KIT



A IMN4164
16 Instrument Infinity Series™ Cassette, Purple
 IMS® cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.



B XP23/86
23/CP-8 Expro
 Double-ended instrument with an explorer on one end and a probe (color-coded) on the other. ergonomic, Satin Steel™ handle.



C H2
2 Kelly Hemostat, Curved
 Used to clamp off blood vessels, remove small root tips and grasp loose objects. Featuring kelly curved beaks.



D PLGGR1X
Grisdale Bone Condenser 8-10-13-15-18, Black Line
 Used during grafting procedures with 8-11-13-15-18 mm markings. Featuring a non-reflective surface, improved contrast between tissue/scissor.



E CM2/4SX
2/4 Molt DE Curette, Black Line
 Designed for curettage, cyst removal and tooth socket debridement with 3.6mm and 7mm working ends.



F 10-130-05
5 Scalpel Handle, Round, Straight
 Used for cutting gingival tissue and making surgical incisions.



G PPR3X
PR-3 Prichard Periosteal, Black Line
 Used for reflecting and retracting the mucoperiosteum. Featuring a large flat blade and narrow sharp blade.



H NH5037
Halsey Perma Sharp® Needle Holder, Serrated
 Used to grasp and guide the needle when suturing. Featuring carbide inserts, finger ring style. Use on all types of sutures.



I S5055
Metzenbaum Perma Sharp® Scissors, Curved/Blunt
 Used for trimming tissue or cutting sutures. 14.5cm



J TP5090
DeBakery Perma Sharp® Tissue Pliers, 6"
 Used to grasp and stabilize soft tissue during suturing.



K DP1
#1 Dressing Pliers, Perry
 Used to grip or transfer materials in and out of the oral cavity. Features a serrated tip to provides a positive grip.



L PMPIX
Membrane Placement Instrument, Black Line
 Micropoint tip for precise placement and manipulation of membranes.



M MCUPE
Measuring Cup 25 mL w/ Marking
 Secure and maintain open position of scissors, hemostats, and other hinged instruments in lid of cassette.



N SYBG
Syringe Bone Grafting 7mm
 Used during grafting procedures. Inner Diameter 5.8 mm. Load from tip to prevent clogging.



O CL856
85 Lucas DE Surgical Curette
 Designed for curettage, cyst removal and tooth socket debridement. 2.5mm



P PJO/16
S13K/TG Jovanovic 4mm Chisel
 For removing and recontouring bone.