



HU-FRIEDY SURGICAL KIT

MICROSURGICAL KIT, BLACK LINE



PART CODE	DESCRIPTION
10-130-70	7 Mini Scalpel Blade Holder
S13K/13KL9	13K/13KL DE Surgical Chisel, #9 Handle
C36/379	36/37 Back-Action DE Chisel, #9 Handle
8-905DD	Precision Dressing Forceps, Diamond Dusted
NH5024CORT	Cortellini Micro Castro Needle Holder
PPBUSERX	Buser Periosteal Elevator, Black Line
SPVX	Castro Microsurgical Scissors, Black Line
IMN4121	12 Instrument Infinity Series™ Cassette, Gray
PCPUNC156	CC SE UNC Probe, 1-15, #6 Handle
MIR5	5 FS CS Mouth Mirror
MH6	Cone Socket Mirror Handle, #6 Handle
IMNMICSURX	Complete Microsurgical Kit, Black Line

MICROSURGICAL KIT, BLACK LINE



* pictured in IMN4121 - 12 Instrument Signature Series® Cassette, Gray

MICROSURGICAL KIT BLACK LINE



- A 10-130-70**
7 Mini Scalpel Blade Holder
Designed to hold micro scalpel blades for minimally invasive procedures and microsurgery



- B S13K/13KL6**
13K/13KL DE Surgical Chisel
Used to recontour periodontal bone or to remove remaining tooth fibers during periodontal surgery with curved working ends 3.2mm and 4.9mm wide.



- C C36/376**
36/37 Back-Action DE Chisel
Designed for removing bone adjacent to the tooth without causing trauma and is especially useful on the distal of last molars with working ends 3.8mm wide.



- D 8-905DD**
Precision Dressing Forceps, Diamond Dusted
Used to grasp and stabilize soft tissue, 7" length.



- E NH5024CORT**
Cortellini Micro Castro Needle Holder
Used to grasp and guide the needle when suturing.



- F PPBUSERX**
Buser Periosteal, #6 Handle
Featuring a sharp round blade and narrow pointed end for interproximal use. For reflecting and retracting the mucoperiosteum.



- G SPVX**
Microsurgical Scissors
Used for trimming tissue or cutting sutures. 16.5cm



- H IMN4121**
12 Instrument Infinity Series™ Cassette, Gray
IMS® cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, and can be specialized to meet your needs.



- I PCPUNC156**
CC SE UNC Probe, 1-15, #6 Handle
Single-ended, color-coded probe. ergonomic, Satin Steel™ handle. University of North Carolina.



- J MIR5**
5 FS CS Mouth Mirror
Single-sided Cone Socket mirror. 15/16 in. (24mm) diameter.



- K MH6**
Cone Socket Mirror Handle, Satin Steel
Used with Single and Double-Sided Cone Socket Mirrors. The blunt end of the mirror handle is ideal for testing sensitivity to percussion and assessing tooth mobility.