

BACK TO PRACTICE

ADDITIONAL 20% OFF AEROSOL MITIGATION KITS

Take advantage of our exclusive offer for the essential aerosol mitigation products you need, bundled for your convenience.

The ADA, CDC, and ADHA recommend prioritizing the use of hand scalers as an alternative to ultrasonic and power methods. Many studies have shown that the same clinical outcomes are achieved with hand scaling as ultrasonics, allowing you to provide enhanced patient care while mitigating aerosol risks. Our EverEdge™ 2.0 Scalers are 72% sharper than the next leading competitor making this essential scaler kit a perfect choice for dental hygiene procedures.



| | |
|---|---------------------|
| Perio Recovery Kit | RECPERIO-KIT |
| Perio Recovery Kit with HVEsolo™ Tips Includes everything from the Perio Recovery Kit, PLUS adds a pack of 50 HVEsolo™ Disposable Evacuation Tips (1808S) | RECPSOLO-KIT |

Your Essential Hand Scaling Kit For Aerosol Mitigation Includes:



H6/H7 Hygienist Scaler | SH6/79E2

Super-thin, contra-angle design. Double-ended, easy-grip handle. Very effective for interproximal areas. Shank angulation provides accessibility to all anterior and premolar surfaces.



204S Sickle Scaler | S204S9E2

For posterior supra-gingival interproximal scaling.



5/6 Barnhart Curette | SBH5/69E2

Barnhart Universal Curette blade designed with two cutting edges and a rounded toe for use on all surfaces.



11/12 Gracey Curette | SG11/129E2

Angulated to reach mesial surfaces of posterior teeth. Area-specific allows for deep scaling, root planing and periodontal debridement. Offset blade provides a perfect working angulation for the tooth surface.



13/14 Gracey Curette | SG13/149E2

Angulated to reach distal surfaces of posterior teeth. Area-specific allows for deep scaling, root planing and periodontal debridement. Offset blade provides a perfect working angulation for the tooth surface.



HVEsolo™ Tips | included only with the Perio Recovery Kit with HVEsolo™ Tips (RECPSOLO-KIT)

Short (2.5" length), lightweight shaft allows for solo practitioner use of high volume evacuation (HVE) in non dominant hand.

Order Now

Contact a **HuFriedyGroup** representative 1-800-Hu-Friedy or your participating distributor.

©2020 Hu-Friedy Mfg. Co., LLC. All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, or Palmero Healthcare., its affiliates or related companies, unless otherwise noted. Offer(s) valid in the 50 United States and District of Columbia. Offer(s) valid June 1 - December 31, 2020. 1. Aerosols and Splatter in Dentistry; Harrell and Molinari; JADA 135; April 2004 HFL-445/0820

Hu-Friedy and Palmero are now proud members of

HuFriedyGroup
The Best In Practice