



**VISTACLEAN™ IRRIGANT
SOLUTION CONCENTRATE**

by Hu-Friedy

IMS-1452

Hu-Friedy Mfg. Co., LLC

3232 N. Rockwell St. | Chicago, IL 60618 | USA

Hu-Friedy

VistaClean is an aqueous cleaner derived from natural citrus botanicals designed for use in dental unit water lines. It may be used continuously or intermittently as desired. When used with **VistaClean®** Dental Water Line Treatment System, follow directions in system manual. Avoid contact with eyes. Keep from freezing or heat. NEVER USE FULL STRENGTH.

VistaClean solution d'irrigation concentrée par Hu-Friedy

VistaClean est un produit nettoyant aqueux dérivé de citrus naturel d'origine végétale, conçu pour être utilisé dans les conduites d'eau des unités dentaires. Il peut être utilisé soit en continu, soit par intermittence. Lorsque ce produit est utilisé avec le système **VistaClean** de traitement des conduites d'eau des unités dentaires, suivez les instructions du manuel du système. Éviter le contact avec les yeux. Ce produit ne doit pas geler; tenir à l'écart de la chaleur. NE JAMAIS UTILISER EN CONCENTRATION INTÉGRALE.

1 Fl. Oz. (29.57ml)

1 oz liquide (29,57 ml)

Best if used before For Vendor
Use Only
Utiliser de préférence avant

 K7036/0413