

Wrap Size Recommendations

	10" x 10"	15" x 15"	20" x 20"	24" x 24"
Signature Series [®] Cassettes				
Exam, Double-Decker™ 8-Instrument, Ultrasonic Insert, Symmetry IQ [®] & Handpiece Cassette	X	X		
Space Saver Seven™ Cassette	X	X		
Double-Decker 14-Instrument Cassette Δ		X	X	
Small 8 or 10 Instrument Cassette Δ		X	X	
Signa-Stat or Double-Decker 20-Instrument Cassette			X	
Large 16 or 20 Instrument Cassette			X	
Oral Surgery Cassette				X
Orthodontic Cassette Box				X
Large Thin Orthodontic Cassette			X	
Small Thin Orthodontic Cassette			X	
Double-Decker Large Ortho Cassette			X	

Δ The Double-Decker 14-instrument cassette and Small cassettes are best with the 15" x 15" sterilization wrap. However, if an office is trying to purchase only one size, the 20" x 20" size is recommended.