

N95 Qualitative Fit Test Kit

For use with any disposable filtering facepiece and negative-pressure half-face respirator¹

QLFT 50 Saccharin N95 Fit Test Kit | 065000-HF



The Occupational Safety and Health Administration (OSHA) requires that fit tests be conducted for those who need to wear respiratory protection, including N95 respirators, to ensure proper function and protection². This Qualitative Fit Test (QLFT) Kit provides the ability to perform pass/fail fit tests that assesses the adequacy of respirator fit based on the individual's response to a test agent.

Features and Benefits

- Simple protocol that is easy to administer
- Includes solution for approximately 55 tests
- OSHA and National Institute for Occupational Safety and Health (NIOSH) compliant

Kit Component	Kit Quantity
A Hood	1
B Saccharin Fit Test Solution, 60mL	1
C Saccharin Sensitivity Solution, 60mL	1
D Nebulizers	2
Replacement Aerosolizing Nozzles	4
Cleaning Tools	2
Instructional USB Flash Drive	1



Refill Component	Refill Part Code	Refill Quantity
Hood	065400-10-HF	10
Saccharin Fit Test Solution, 60mL	065200-HF	6
Saccharin Sensitivity Solution, 60mL	065100-HF	6
Nebulizers	065300-HF	3
Bitrex® Sensitivity Solution, 60mL	066100-HF	6
Bitrex® Fit Test Solution, 60mL	066200-HF	6

For more information call 1-800-Hu-Friedy, contact your authorized HuFriedyGroup representative, or visit Hu-Friedy.com

¹) Respirator not included in Fit Test Kits ²) OSHA Standard 1910.134 App A

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC its affiliates or related companies, unless otherwise noted. Bitrex is the registered trademark of Johnson Matthey/MacFarlan Smith.

©2020 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-481/0820