

TECHNICAL DATA SHEET



PRODUCT

Crosstex CSR Sterilization Wrap

PRODUCT OVERVIEW

- Can be used in all types of sterilization centers and surgical unit packaging.
- Validated for all steam sterilization; dynamic air removal (pre vac and steam flush pressure pulse); and gravity displacement equipment and EtO; gamma irradiation and E-beam sterilization technologies
- Can be used for inner and outer wrap (double layer) for a varying cassette sizes
- CSR Wraps are single use item

MATERIALS

- 52 gsm cellulose fibers reinforced with synthetic binders
- Cellulose pulp >80%*
- Binders <20%
- Sizing Agent 0.3%
- Wet Strength Agent 0.3%

Note: (*) Petroleum Free – Manufactured from environmentally friendly cellulose material made from wood pulp.

DETAILS

- Excellent bacterial barrier to both airborne and water-borne contamination
- Outstanding liquid repellency – water, alcohol, iodine
- Exceptional tensile strength – wet or dry
- Even texture to allow immediate visual detection of small pinholes or tears
- Meets FFA (Flammable Fabrics Act) requirements
- Complies with EN ISO 11607-Part 1 and EN 868 – Part 2 All validation

**PRODUCT/
PACKAGING**

Product is available in the following configurations

Product Code	Description
SW12	12" x 12" – 1,000/Case
SW15	15" x 15" – 500/Case
SW18	18" x 18" – 500/Case
SW20	20" x 20" – 500/Case
SW24	24" x 24" – 500/Case
SW30	30" x 30" – 250/Case

Product Code	Description
SW36	36" x 36" – 250/Case
SW40	40" x 40" – 250/Case
SW45	45" x 45" – 250/Case
SW48	48" x 48" – 100/Case
SW54	54" x 54" – 100/Case

STORAGE REQUIREMENTS

Store in a cool, dry location always from strong sunlight or UV light (30 to 60% RH; 50°F to 104°F / 10°C to 40°C).

SHELF LIFE

When properly sealed, sterility is maintained up to (5) five years provided the integrity of the wrap has not been compromised.

0818 DTDS00327 (MKT769)

All product names are trademarks of Crosstex International, Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted.