

Premium™ Saliva Ejectors

Advantage Saliva Ejectors

INSTRUCTIONS FOR USE

1. Attach ejector head to hand piece.
2. Place in mouth and bend to desired position.
3. After use, dispose of the ejector.



Ref.#'s ZCB, ZCW, ZWW, ZPP,
ZCMC, ZBBI, ZCBI, ZCCI, ZWWI



Crosstex International, Inc.
6789 W. Henrietta Road, Rush, NY 14543 USA
585.359.0130 | 631.582.6777 | crosstex.com



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Not made with natural rubber latex.

Made in USA
Z-IFU-EN / REV B - 03/2022

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.
©2022 Crosstex International, Inc. All rights reserved.

CROSSTEX